



SprintPack™ Lithium-Ion  
Power System

# Unlimited, uninterrupted hot-swappable power

The SprintPack Lithium-Ion Power System is equipped with advanced technology designed to empower its users with comfort, confidence and convenience.

Lithium-ion technology provides the greatest degree of power density in the lightest-weight package of any battery type. This provides LTV® ventilator users longer battery power in a smaller and lighter package. The SprintPack Lithium-Ion Power System also provides a higher number of full recharging cycles that extends the life of the battery to provide greater long-term value.

## Unlimited, uninterrupted hot-swappable power

When using the SprintPack Lithium-Ion Power System, you may take out one battery at a time, leaving the other charged battery in place without interrupting power. Upon replacing a spent battery with a charged power cell, the SprintPack power manager automatically switches to the remaining, installed charged power cell. This allows you to swap spent cells with recharged cells almost indefinitely without losing power to the ventilator.

## Uninterruptable power supply functionality

The SprintPack Lithium-Ion Power System operates as an uninterruptable power supply for the LTV ventilator and works in line with the ventilator and power supply, simultaneously recharging and powering the ventilator. If the power from the wall outlet or other source is interrupted, the SprintPack Lithium-Ion Power System immediately takes over to power the ventilator without interruption or any action by the user or caregiver.

## Lighted power gauge

SprintPack power cells use a lighted gauge to indicate the amount of charge remaining in the cell. SprintPack batteries can be charged from AC or DC power sources, such as in an automobile, in an emergency service vehicle or when connected to a wheelchair battery. The AC power adapter can recharge from wall outlets in line with the ventilator or separately to prepare for the next use.

## Key features include:

- Hot-swappable
- Multiple recharging options
- Uninterruptable power supply
- LED charge level indicators
- Lightweight lithium-ion technology



## Physical specifications

<b>Weight</b>	4.5 lb (2.04 kg)
<b>Height</b>	1.65" (4.19 cm)
<b>Width</b>	6.70" (17.02 cm)
<b>Depth</b>	10.85" (27.56 cm)

## Power specifications

<b>Input voltage</b>	11 to 29 VDC
<b>Output voltage</b>	11 to 15 VDC
<b>Operating time*</b>	2.5 hr, 1 battery 5 hr, 2 batteries
<b>Recharge time</b>	3 hr, 1 battery 5 hr, 2 batteries
<b>Battery chemistry</b>	Lithium-ion
<b>Cycles</b>	400 full cycles at 70 °F (21 °C)



## Environmental specifications

<b>Operating temperature</b>	41 to 104 °F (5 to 40 °C)
<b>Storage temperature</b>	-4 to 122 °F (-20 to 50 °C)
<b>Operating humidity</b>	15 to 95% relative
<b>Storage humidity</b>	10 to 95% relative

## Power indicators

<b>Charge status</b>	5 segment LEDs
<b>10% remaining charge</b>	1st LED flashes yellow at 2 Hz
<b>Recharging</b>	LED blinks current capacity at 1 Hz
<b>Fully recharged</b>	All 5 LEDs illuminated

## Standards compliance

<b>Shock and vibration</b>	MIL-STD-810F, IEC 68-2-27, IEC 68-2-34, IEC 68-2-6
<b>Water ingress</b>	IEC 601-1, Clause 44.3
<b>EMC compatibility</b>	IEC 60601-1-2 Emissions, immunity and ESD
<b>RF emissions</b>	RTCA DO-160D Section 21.4 Cat. M
<b>ESD protection</b>	ASTM F1246-91 4.14.3.1
<b>Safety</b>	IEC/UL 60601-1
<b>Shipping/transport</b>	Normal cargo (< 8 g lithium per pack)

\* At nominal ventilator settings



For more information, please contact Customer Service within the United States at: **800.754.1914** and International at: **+1 763.398.8300**.

#### GLOBAL HEADQUARTERS

Vyaire Medical, Inc.  
26125 North Riverwoods Blvd  
Mettawa, IL 60045  
USA



Vyaire Medical, Inc.  
26125 North Riverwoods Blvd  
Mettawa, IL 60045  
USA



EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

#### AUSTRALIAN SPONSOR

Vyaire Medical Pty Ltd  
Suite 5.03, Building C  
11 Talavera Road  
Macquarie Park, NSW 2113  
Australia



[vyaire.com](http://vyaire.com)

**For global distribution.**

Trademarks are the property of their Respective owners.

© 2020 Vyaire. Vyaire, the Vyaire logo and Vyaire SprintPack Lithium-Ion Power System are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-2000161