

AirLife™ brand Misty Max 10™ nebulizer

Respiratory care products and services

David A. Rivera, Research and Development

December 18, 2002

Technical data

The purpose of this document is to provide the in vitro performance characteristics of the AirLife™ brand Misty Max 10™ nebulizer. In vitro measurements will help guide the end user as to the likelihood of effective delivery of inhaled medication to the lower respiratory tract.¹

Particle size

Methods

An Andersen eight-stage cascade impactor (ACI) with USP inlet² was used and assayed with a spectrophotometer (224 nm). Particle size characterization was performed using 3 mL of 0.083% albuterol sulfate in normal saline. Cascade impaction was chosen for the following reasons:

1. Cascade impaction measures aerodynamic diameter directly, which accounts for the density and irregular shape of drug particles. It is believed that aerodynamic diameter more accurately predicts the behavior of aerosol as it is delivered into the patient's lungs.³
2. There is more historical data on particle size measurement using cascade data than any other method. Relative comparison with historical data can be readily made.
3. Cascade impaction is one of the USP methods for characterization of particle distributions.²

Performance characteristics	7 L/min flow	8 L/min flow	10 L/min flow
MMAD Mass median aerodynamic diameter (µm)	1.66 ±0.06	1.48 ±0.02	1.35 ±0.04
GSD Geometric standard deviation	1.45 ±0.04	1.44 ±0.02	1.45 ±0.03
Respirable fraction % < 4.7µm (by mass)	85.2 ±4	91.4 ±5	93.9 ±6
Output rate (g/min)	0.313 ±0.022	0.312 ±0.017	0.338 ±0.014
Treatment time 3 mL albuterol solution (min)	8.66 ±0.60	8.77 ±0.55	8.24 ±0.71
Residual (g)	0.84 ±0.06	0.84 ±0.07	0.84 ±0.07
Operating pressure (psi)	20.8 ±0.2	25.8 ±0.3	36.0 ±0.0

Data presented: Mean ±standard deviation.

Sufficient sample sizes were tested to obtain a 95% confidence interval on the mean less than or equal to: ±9% respirable fraction and ±4% for all other characteristics.

Test conditions

Particle size testing was conducted with components included with AirLife brand Misty Max 10 nebulizer 002438: T-adapter, mouthpiece, 7' oxygen tubing and 6" flex tube. The nebulizer was driven by compressed air (regulated to 50 psi) at 7, 8 and 10 L/min. flow. Ambient air, which flows through the open end of the T-adapter, was controlled at 50% ±10% relative humidity.

Definition of parameters measured

1. MMAD: Mass median aerodynamic diameter

This is a measure of central tendency of the size distribution of aerosol particles. It is the diameter, in micrometers, of which 50% of the mass of aerosol is larger and 50% is smaller.

2. GSD: Geometric standard deviation

This is a measure of the width of the size distribution of aerosol. For log normal distributions, it is calculated as follows:² $GSD = (84.13\% \text{ diameter} / 15.87\% \text{ diameter})^{1.2}$

Respirable fraction

Respirable fraction is the percent of aerosol generated by mass that falls below 5 μ MMAD. It has been reported that particles less than 5 μ will penetrate beyond the upper airways and deposit into the tracheobronchial and pulmonary regions of the lung.⁴ This measurement is often used to describe the quality of aerosol.⁵ The closest cascade plate cutoff point of 4.7 μ was utilized to quantify the mass below 5 μ .

Output rate

Methods

Nebulizers were weighed before and after aerosolizing a fill volume of 5 mL normal saline for 3 min. Performance was reported in the average weight loss per minute (g/min).

Test conditions

Testing was done with nebulizers randomly pulled from production lots. The nebulizer was operated with compressed air, and normal saline was the nebulized solution.

Residual and treatment time

Methods

The measurements of residual and treatment time will depend on the criteria used to end the treatment. In clinical practice, it is common to tap the nebulizer near the end of the treatment to maximize the amount of drug nebulized. To simulate this practice, the nebulizer was periodically tapped until tapping resulted in an aerosol production of less than 5 sec.

Residual was calculated by subtracting the weight of the nebulizer after the trial from the weight of the nebulizer empty. Treatment time was the total time required to nebulize 3 mL of solution.

Test conditions

Residual and treatment time testing was done with nebulizers randomly pulled from production lots. The nebulizer was operated with compressed air, and normal saline was the nebulized solution.

Operating pressure

Methods


Pressure was obtained from a T-connection between the compressed air flow meter and the nebulizer.

Test conditions

Testing was done with nebulizers randomly pulled from production lots. Pressures at 7, 8 and 10 L/min. flow were evaluated while nebulizing 3 mL of normal saline.

References

- 1 Consensus Statement: Aerosols and Delivery Devices. *Respiratory Care*, 2000, 45(6).
- 2 Aerosols. U.S. Pharmacopia, 23(601).
- 3 Hess, D. Medication nebulizer performance. *Chest*, 1996, 110(2): 498-505.
- 4 Laube, B. In vivo measurements of aerosol dose and distribution: Clinical relevance. *Journal of Aerosol Medicine*, 1996, 9(1).
- 5 Dolovich, M. Influences of inspiratory flow rate, particle size and airway caliber on aerosolized drug delivery to the lung. *Respiratory Care*, 2000, 45(6).

 CareFusion
22745 Savi Ranch Parkway
Yorba Linda, CA 92887
800.231.2466 toll-free
714.283.2228 tel
714.283.8493 fax

 CareFusion France 205 S.A.S
5 rue Charles Lindbergh
F-44110 Châteaubriant, France
+33 240 813 804 tel
+33 0 240 815 634 fax

© 2011 CareFusion Corporation or one of its subsidiaries. All rights reserved. AirLife and Misty Max 10 are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. RC2342 (0211/1000) L3370 Rev. A

CareFusion
Yorba Linda, CA

carefusion.com

