


AioCare Spirometry Better Diagnostics Better Lives


A completely mobile and flexible approach



FLOW/VOLUME



VOLUME/TIME



AioCare system specifications

Measurement



SVC Slow spirometry, FVC - Forced spirometry, BRT - Bronchodilator Responsiveness Test, SpO2 & Heart rate, Peakflow Diary



Parameters		
Slow parameters	VC, IC, ERV,	VT, IRV
Forced parameters	FEV1, FVC, FEV1/FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEV6, FIVC, PIF, FIF25, FIF50, FIF75	
Technical parameters	VPTEF	Volume to Peak Tidal Expiratory Flow
	ERV	Expiratory Reserve Volume
	TPTEF	Time to reach Peak Tidal Expiratory Flow
	TPTEF/TE	Time to reach Peak Tidal Expiratory Flow as a proportion of Total Expiratory time
	TPEF	Time to Peak Expiratory Flow. Time from the start of the forced exhalation to the point of Peak Expiratory Flow
	RT	Rise Time. Time required for a signal to change from 10% to 90% of the Peak Expiratory Flow
	TPTEF VC	Ratio of Volume to peak expiratory flow to total expiratory volume
	VPTEF/VE	Ratio of Volume to peak expiratory flow to total expiratory volume
	VT	Volume Tidal
	FET	Forced Expiratory time
	BEV	Back Extrapolated Volume

Flow measurement	
Sensor type for flow measurement	Thermal
Spirometric flow measurement range	0-16 l/s
Flow accuracy	±5% or 200 mL/s
Resistance	<0.5 cm H2O/L/s
Volume range	0-8 litres
Volume accuracy	±2.5% or 50 ml, whichever is greater
Linearity	2.5%


Volume integration		_
Flow measurement resolution	Measured 5 ml/sec, usable 10 ml/sec	/
Accuracy/Repeatability	Standard: ATS/ERS 2019	
Automatic BTPS correction	Built-in sensors for measuring temperature, pressure & humidity	
Determination of t0	Algorithmic	
Expiratory impedance	<0,15 kPa/(I/s) at 14I/s	

Technical		
Protection of the casing against water ingress, according to IEC 60529 (spirometer elements)	IP 22	
Communication	Bluetooth 4.0. Low Energy	
Bluetooth frequency	2.4-2.48 GHz	
Measurement frequency	100 Hz	
Internal power supply	Battery (LiPo 3.7 V)	
50 mA power consumption	50 mA	
Dimensions	118 x 38 x 48 mm	

Weight

0,3 kg

Standars, directives and material o	clearances)
Standards	ATS/ERS 2019, EN 60601-1, EN 60601-1-2, EN 62304, EN 62366, EN ISO 14971, EN ISO 10993-1	
Directives	93/42/EEC amended by 2007/47/EC, RoHS 2011/65/EU compliant	
Market clearances	CE 2294	









No specific requirements Internet connection

Requirements

AIoCare Panel





Highest quality standards



Patient Validation Tests



Multicenter Observation Study



FVC Accuracy



National obstruction Diagnostic Program



Implemented 2019 guidelines



Accuracy Validation Report, May 2020

AioCare + MicroGard™ II Filter Protects patients

High level validated cleaning methods

Protect patients, staff and device from cross-contamination by using the MicroGard™ II in-line filter.

MicroGard II filters provide 99.999% viral and bacterial efficiency against cross contamination (Nelson Test Report 1003754).



The AioCare spirometer and antibacterial/viral filter MicroGard II respectivley manufactured by HealthUp and Vyaire were tested according to ATS/ERS Standardization of Spirometry 2019.

This testing is to verify the quality of the results considering accuracy, repeatability, linearity and resistance to flow of the combination.

All 50 AioCare spirometers (which have been tested using MicroGard II antibacterial/viral filters) have met the full criteria described in ISO23747:2015 and ISO26782:2009 standards. The variety of waveforms in both standards encompass the characteristics seen in the population of patients.

- Accuracy, repeatability and linearity for waveforms C1-C11 (applies to 26782:2009 standard) are within the permissible error range.
- Accuracy, repeatability and linearity for waveforms C12-C13 tested on heated and humidified air as well as impedance test (applies to 26782:2009 standards) are within the permissible error range.
- Accuracy, repeatability, linearity for profiles A and frequency response (applies to 23747:2015 standard) are within the permissible error range.
- Accuracy, repeatability, linearity for profiles A 300I/min and 600I/min as well as impedance tests (applies to 23747:2015 standard) are within the permissible error range.

ISO 26782:2009

Anaesthetic and respiratory equipment — Spirometers intended for the measurement of time forced expired volumes in humans

ISO 23747:2015


Anaesthetic and respiratory equipment — Peak expiratory flow meters for the assessment of pulmonary function in spontaneously breathing humans

Server / Backend



AioCare infrastructure is supported by Amazon AWS platform. AWS provides physical security and we are operating from Central Europe region that has its DataCenter located in Frankfurt, Germany. AWS infrastructure meets all the required standards: ISO/IEC, CSA/CCM, ITAR, CJIS, HIPAA and IRS. The code runs on AWS Elastic Container Service infrastructure and its access and management is fully automated using PK/SSH security. Any manual intervention is only available to AWS System Administrator. In addition, access to any Data is restricted for specific IP addresses. Only services necessary for user operation are available from the internet, using always HTTPS protocol. The many modules that compose AioCare Solution are written in PHP v8.1 and node.js v18.1, using a Relational DataBase (MySQL8.2).

Access to data through the API is possible only after providing correct user credentials, that will grant the user a JWT token, that is only valid for the next 2hours on Web and 24hours on Mobile, after that period the user needs to authenticate again.



The app allows 3 different roles, and requires different layers of authentication:



System Administrator

- Email/Password
- 2FA using email code
- IP restriction



Doctor

- Email/Password
- 2FA using email code



Patient

- Email/Password

All authentication data (ie: password) is stored encrypted and its access is fully restricted and audited. Any manually required intervention needs to be authorized by the DB administration and is fully audited. Data backups are performed daily and are automated and stored for the previous 7 days.

REFERENCES

1. Based on the Bio Burden DIN EN ISO 11737-1: Report 18AA0088

The contents of this publication may differ from the current approval of the product or service in your country.

GLOBAL HEADQUARTERS

Vyaire Medical, Inc. 26125 N. Riverwoods Blvd. Mettawa, IL 60045 USA

AIo**C**are

HealthUp S. A.
Twarda 18
00-375 Warsaw
Poland

C€ 2274

MICROGARD™ II FILTER SENTRYSUITE™ SOFTWARE

Vyaire Medical GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany

C € ₀₁₂₃

For international use only.

Trademarks are the property of their respective owners.

© 2023 Vyaire Medical, Inc. or one of its affiliates. All rights reserved. Vyaire and the Vyaire Medical logo are trademarks or registered trademarks of Vyaire Medical, Inc. or one of its affiliates. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-INT-20000237 | 2.0

