

Product Security Bulletin

IMPORTANT INFORMATION – PLEASE READ AND KEEP

BlueKeep Security Vulnerabilities

Date: 2019-05-21

DocumentID: SEC-04

Background

A security vulnerability in the Remote Desktop Protocol (RDP) known as **BlueKeep** ([CVE-2019-0708](#)) has received significant news coverage for its potential impact.

What BlueKeep Does

A remote code execution vulnerability exists in Remote Desktop Services – formerly known as Terminal Services – when an unauthenticated attacker connects to the target system using RDP and sends specially crafted requests. This vulnerability is pre-authentication and requires no user interaction. An attacker who successfully exploited this vulnerability could execute arbitrary code on the target system. An attacker could then install programs; view, change, or delete data; or create new accounts with full user rights.

What Operating Systems Are Affected

Affected operating systems include Microsoft Windows XP, Windows 7 SP1, Windows Server 2003, Windows Server 2008 and Windows Server 2008 R2.

Windows 8 and Windows 10 are **not** affected by BlueKeep.

Details regarding the vulnerability and the Microsoft's impacted products and responses may be found [here](#). Proof of Concept exploit code has been authored but is not yet publicly available.

Product Security Bulletin

IMPORTANT INFORMATION – PLEASE READ AND KEEP

Response

How Microsoft Is Responding

Microsoft has [communicated](#) the issue to its customers on May 14th 2019. Microsoft strongly recommends that all affected systems are patched as soon as possible. [Patches](#) have been made available to download – also for [Windows XP and Windows Server 2003](#) which no longer receive mainstream support.

Affected Products

Vyaire's medical devices affected include all Windows-based Respiratory products which use an affected operating system.

Vyaire has tested all affected products for compatibility with the patches and has found **no** issues. Vyaire recommends to install the patches as soon as possible.

Mitigations & Compensating Controls

Ensure that the patches made available by Microsoft have been applied to the operating system of your computers:

- Windows 7, Windows Server 2008: [Microsoft Advisory CVE-2019-0708](#), published May 14th 2019
- Windows XP, Windows Server 2003: [Microsoft Customer guidance for CVE-2019-0708](#), published May 14th 2019

Generic controls

- Ensure your data has been backed up and stored according to your individual process and that your disaster recovery procedures are in place.
- Update your anti-virus and malware protection, where available.

For product or site-specific concerns, contact your Vyaire service representative.

For more information on Vyaire's proactive approach to product security and vulnerability management, contact us under: <http://www.vyaire.com/productsecurity>