


# SentrySuite<sup>™</sup>

software platform for intelligent respiratory diagnostics



# SentrySuite™ 360° – one screen, all in




# SentrySuite™ supports a comprehensive list of PFT and CPET measurements



Whether the **MasterScreen<sup>™</sup>** series or the new **Vyntus<sup>™</sup>** series, they are all **powered by SentrySuite** – from small PC spirometers to high-end PFT to ECG and CPET devices.

All PFT measurement programs show the same graph and table arrangement with a consistent user interface. Intuitive operation and quick learning of new measurements are guaranteed.





# Measurement capabilities:

**SB Diffusion Real-time** 


Slow Spirometry	SB Diffusion Intrabreath	Compliance
Forced Spirometry	NO Diffusion	Resting ECG
MVV	Rocc	Stress ECG
Bronchial Challenge Testing	Impulse Oscillometry (IOS)	CPET breath-by-breath
Body Plethysmography	Rhinomanometry	<b>CPET</b> mixing chamber
FRC N <sub>2</sub> Washout	P0.1	CPET High/Low FIO <sub>2</sub>
FRC HE Rebreathing	MIP/MEP	REE Canopy


**SNIP** 

# **Guidance and Quality**

Coaching for operator and patient

**Check marks** for spirometry end-of-forced-expiration criteria





**Breath-hold timer** and mouth pressure control

> Online textual guidance ideal for new personnel

## **Animation incentives**


Choose from 10 different patient animations for forced spirometry.












Adjust the target scale to your patient needs.

## SentrySuite<sup>™</sup> Quality Assurance-approved results

SentrySuite offers many tools to approve the quality and accuracy of your measurements.

> With SentrySuite Quality Management you can validate your measurement device using calibrated syringes, FRCpleth or DLCO simulators. Pre-defined special SeS QM patients automatically switch-off BTPS compensation and correct appropriately for mechanical dead space.




SentrySuite supports biologic control testing with a separate **BIO QC database**. The BIO QC data is easily identified and separated out for trending and further analysis.

According to the ATS/ERS guidelines, tests should be performed on a "standard subject".

After each trial, SentrySuite provides feedback on the acceptability and repeatability of the measurement, checking the results against the ATS/ERS guidelines or other defined standards.

- ATS/ERS Spirometry Guidelines 2019
- ATS/ERS Diffusion Guidelines 2017
- ATS/ERS Reporting Guidelines 2017












# **Networking**


The united data of SentrySuite™

## **Central Patient Database and HIS Interface**

Server - SQL unlimited - HL7 - web-based - HIS - DICOM - GDT



# Save time with the workflow driven physician's review tool




# With SentrySuite™, technician and physician work hand in hand

The technician performs the measurement and saves the proper report.

The physician sees all visits with their associated reports for interpretation.

Flexible workflow status design allows for optional steps such as a preliminary interpretation by a fellow, with subsequent confirmation or revision by the attending physician.

Apply filters to create individual work lists.



## **Review Tool**

A standard with SentrySuite™

- Ultra-fast report view with results on the left and interpretation on the right
- Customized interpretation templates and macros
- Electronic signature
- Automated Send-To-HIS functionality
- Status management
- Also available as web-based solution named SentrySuite Mobile Review


## **Artificial Intelligence meets PFT interpretation**

Seamless workflow integration

**ArtiQ.PFT** fully automates the interpretation of pulmonary function tests and provides Al-based disease probabilities to support diagnosis. Allowing pulmonologists to focus on clinical decision making and patient consultation.



- spend your time wisley
- boost consistency
- maximize your diagnostic accuracy
- define your next steps




# Tools go hand in hand with the workflow

		Pred Auth	Pred	Pred LL	Best	%(Best/Pred)	Tr 1	Tr 2	Tr 3	Z-Score	5 -4 -3	Z-Score
FVC	L	Quanjer GLI (2012)	5.06	4.05	4.84	96	4.84	4.78	4.68	-0.37		8
FEV 1	L	Quanjer GLI (2012)	4.13	3.28	4.07	98	4.07	4.02	3.97	-0.12		6
FEV 1 % FVC	%	Quanjer GLI (2012)	81.86	71.28	84.15		84.15	84.10	84.82	0.40		8
PEF	L/s	ECCS (1993)	9.33	7.34	11.09	119	11.09	10.87	10.53	1.46		
MFEF 75/25	L/s	Quanjer GLI (2012)	4.14	2.48	4.26	103	4.26	4.23	4.18	0.10		0
Quality Grade FVC ATS 2019					Α							
Quality Crade EEVA ATC 2040					A							

## Scoring

SentrySuite offers all important parameter related scores; %Predicted, Lower Limit, Upper Limit, Z-Score, Standard Deviation




## Classification

Each measurement program shows a classification bar indicating the possible limitation



## **Interpretation Templates**

Create your own templates with macros



## **Automatic Interpretation**

Several authors come standard: Ellis, Jaeger (1994), ATS/ERS 2005, Intermountain Thoracic Society (1984), Methodist Hospital (1994),


IOS, CPET Eschenbacher/Mannina (1990), Fukuda Sangyo



## **Individual Bar Graphs**

Create your own interpretation bar with the flexible Report Designer

# Trending-monitor patient results over time




nt		Spiro		
		FVC	FEV1	MMEF
Time	Medication	(L)	(L)	(L/s)
16:26:28		4.85	4.07	3.81
16:32:53	Salbutamol	4.87	4.05	3.70
13:22:44		4.84	4.07	3.99
12:57:24		4.72	4.02	4.26
11:12:17		4.79	4.03	4.19
10:55:39		4.57	3.73	3.64
		4.77	3.99	3.93
		0.10	0.12	0.23
]		2.13	2.94	5.96
	16:26:28 16:32:53 13:22:44 12:57:24 11:12:17	Time Medication  16:26:28 16:32:53 Salbutamol 13:22:44 12:57:24 11:12:17 10:55:39	Time Medication (L)  16:26:28	Time         Medication         (L)         (L)         (L)           16:26:28         4.85         4.07           16:32:53         Salbutamol         4.87         4.05           13:22:44         4.84         4.07           12:57:24         4.72         4.02           11:12:17         4.79         4.03           10:55:39         4.57         3.73           4.77         3.99           0.10         0.12



Trending reports are standard in SentrySuite. Show your long term data with one click for any parameter.

# Respiratory data transfer

Connect your Fenom Pro® and stand alone Vyaire spirometry devices with SentrySuite™.



# Blood gas online connection


Connect your blood gas analyzer with SentrySuite

- Automate import and storing of blood gas measurements into the SentrySuite patient data record
- PFT / CPET and blood gas results on one report ready for interpretation
- Supports many types of blood gas analyzers

# **Questionnaires**

Digitally integrated into SentrySuite $^{\text{\tiny{TM}}}$  – paperless and mobile solutions


Collect patient answers electronically and integrate them seamlessly into the SentrySuite database. The completed questionnaires are saved into the patient's database and can become part of your report. Choose from our standard questionnaires or create your own with the optional Questionnaire Designer with six different question types.




# Offline data input

Offline data input – a standard in SentrySuite

Extend and feed your SentrySuite database with further external measurements like exhaled NO, lactate values and blood gases.



on one common report








Define and deploy tailored research databases, based on SentrySuite patient demographics and visit data, standard questionnaires, predicted information and measurement parameters. Easy **Query Builder** user interface

- Direct output into Microsoft Excel or others
- Profit from various table configurations
- Customize, save and load your queries



SentrySuite provides access to your patient data for clinical research.

From SentrySuite Query Builder to Excel Outcome

	A	В	C	D	E	F	G	H	1	GT	HC	HI	JK	JN	LA	LD	П	LY
	Visit Date	Patient ID	Gender	Race	Test Date	Age	Height	Weight	LevelType	FEV 0.5	FEV 1	FEV 1 % FVC	FVC	FVC in	MEF 25	MEF 50	MEF 75	MFEF 75/25
	27.11.2014	GAB10031963	male	Caucasian	27.11.2014	51	1,74	90	1	3,492	4,255	83,613	5,089	5,484	1,137	4,788	11,253	4,643
	13.04.2013	CFNSPJO750209	male	Caucasian	13.04.2013	38	1,74	80	1	3,129	3,986	84,485	4,718	4,691	1,974	4,819	8,318	4,179
5	27.04.2012	CFNSPJO750209	male	Caucasian	27.04.2012	37	1,74	80	1	3,24	4,139	84,014	4,965	5,002	1,868	4,898	9,183	4,188
7	09.09.2011	CFNSPJO750209	male	Caucasian	12.09.2011	36	1,74	80	1	3,145	4,024	85,226	4,721	4,766	2,015	4,785	8,101	4,339
3	25.01.2011	CFNSPJO750209	male	Caucasian	25.01.2011	36	1,74	80	5	2,984	3,779	91,57	4,127	3,408	2,506	5,183	7,748	4,733
9	25.01.2011	CFNSPJO750209	male	Caucasian	25.01.2011	36	1,74	80	5	3,022	3,873	91,855	4,216	3,622	2,642	5,196	7,761	4,748
	25.01.2011	CFNSPJO750209	male	Caucasian	25.01.2011	36	1,74	80	5	1,424	1,848	91,439	2,021	1,399	1,2	2,424	3,628	2,276
			male				1,74	80	6		3,911	82,642	4,732	4,485	1,74	4,502		



# SentrySuite™ Cyber-security

SentrySuite software protects you against cybersecurity threats by implementing a whole set of technical safeguards.


Data security and access control are the keys to protecting health information. SentrySuite offers several security features and comprehensive user administration.

Three steps to secure your patient data:



## 1. AUTHENTICATION

SentrySuite offers **comprehensive settings for your password protection** including the use of Single Sign On and Active Directory linking through **LDAP**.








## 2. AUTHORIZATION

The SentrySuite program **User Administration** allows easy to **create user and user groups** with their **assigned rights.** Make sure that only authorized users have access to the different sub programs.





#### 3. AUDIT

The program **Log File Viewer** reports all proceeded steps within the SentrySuite Software with their appropriate user.

Date Time User name		Workstation ID	Site	Event ID	Event text	Application	
✓ 2017-02-14 10:59:26 ,070	Dr. Smith	LWDEU073XL56R1	n/a	4021	Startup program	Bodyplethysmography	
√ 2017-02-14 10:59:25 ,250	n/a	LWDEU073XL56R1	n/a	4013	Dr. Smith	Session Manager	
√ 2017-02-14 10:59:24 ,090	Dr. Smith	LWDEU073XL56R1	n/a	4021	Exit program	Spirometry	
2017-02-14 10:20:25 ,880	Dr. Smith	LWDEU073XL56R1	n/a	4021	Startup program	Spirometry	
2017-02-14 10:20:25 ,150	n/a	LWDEU073XL56R1	n/a	4013	Dr. Smith	Session Manager	
2017-02-14 10:20:23 ,800	n/a	LWDEU073XL56R1	n/a	4021	Starting Application 'Spirometry'	SentryDesktop	
2017-02-14 10:20:14 ,030	n/a	LWDEU073XL56R1	n/a	1000	Run Automatic: Completed Task 'Start up database'	SentryDesktop	

#### **GLOBAL HEADQUARTERS**

Vyaire Medical, Inc. 26125 N. Riverwoods Blvd. Mettawa, IL 60045 USA

## ARTIQ.PFT

ArtiQ NV

Boskouter 15

3010 Kessel-Lo

Belgium

 $\epsilon$ 

#### SENTRYSUITE" SOFTWARE

Vyaire Medical GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany

**C** € <sub>0123</sub>

#### FENOM PRO®

CAIRE Diagnosics Inc.
7020 Koll Center Pkwy,
Suite 110
Pleasanton, CA 94566-3107
USA



EC REP

QAdvis EAR AB Ideon Science Park Scheelevagen 17 SE-223 Lund, Sweden

#### For International distribution only.

© 2022 Vyaire Medical, Inc. or one of its affiliates. All rights reserved. Vyaire, the Vyaire Medical logo and all other trademarks are trademarks or registered trademarks of Vyaire Medical, Inc. or one of its affiliates. Vyaire's Medical devices are class I & Ila according to Medical Devices Directive 93/42/EEC or Medical Device Regulation EU 2017/745 as indicated on each declaration of conformity. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labeling. VYR-INTL-1900022 | 2.0

