



Vyntus™ CPX

Cardiopulmonary exercise testing at its finest

VYNTUS™ CPX

 **vyairé**™
MEDICAL

Vyntus™ CPX metabolic cart

Striking the perfect balance between high tech flexibility that is easy-to-learn and easy-to-perform.


The versatile Vyntus™ CPX Metabolic Cart combines Vyair's pioneering experience with the latest innovations, to deliver the newest generation of clinical cardiopulmonary exercise testing. Utilizing highly accurate sensors to collect full breath-by-breath gas measurement data, the system provides helpful guidance and tools to make it easier for technicians. And, the fully workflow-driven evaluation helps simplify and standardize data reporting for physicians.

- **Digital Volume Transducer (DVT)** compact, lightweight design with little dead space and minimal resistance to airflow in both cleanable and single-use, disposable options
- **Automatic volume calibration** ensures consistency, saves time and hassle
- **On-board pulse oximetry** with finger, ear-clip and forehead sensors
- **Powered by SentrySuite™** with cues and guidance during measurement and post-test workflow to help standardize evaluations and reduce time to results
- **Smart tools automate** processes for clinicians such as automatic slope calculations and automatic trending of patient data
- **Utilities for customization** including our extensive, global library of adult and pediatric predicted equations and our comprehensive report generation capabilities



Our CPX module

Packed with innovative, robust technologies and features that simplify and improve both testing and maintenance.



No syringe required! CPX automated volume calibration eliminates time-consuming and technique-dependent syringe calibration.

Interface your ECG for comprehensive CPET testing. Seamlessly.

Our flexible Vyntus™ CPX system integrates several commercially available ECG devices; including GE Healthcare CASE™ Exercise Testing System, CardioSoft™, Mortara, Welch Allyn, PBI and others. Vyntus CPX changes everything by changing nothing on your end.

Introducing Vyntus CPX Big Cinema

Vyntus CPX combined with GE Healthcare CardioSoft ECG creates an all-in-one device for a simplified, space saving solution.

CardioSoft™ ECG

- CardioSoft data automatically transfers to Vyntus CPX
- Acquisition module (CAM 14) connects to a standard PC USB port.
- Full Disclosure stores unfiltered, continuous ECG signals.



Vyntus CPX Standard Display

Whether it's the CASE system from GE Healthcare or other existing ECG systems, Vyntus CPX is the perfect compliment. Our standard single monitor display option allows you to view your gas exchange data while leveraging your already existing monitor for ECG.




Vyntus™ CPX Big Cinema.

All gas exchange and ECG information on a single screen. Just look.

- 1 Display your selected metabolic parameters
- 2 View ongoing performance relative to predicted max values
- 3 Quick patient data access
- 4 Tabs to quickly switch to view different graphics
- 5 Manual override of bike or treadmill protocol
- 6 Real-time data with color-coded maximum predicted ranges
- 7 Countdown to upcoming submeasurement programs showing when next programmed events will occur
- 8 60 second view of patient breathing
- 9 Pertinent ECG data is shared with SentrySuite automatically



9 panel graph is viewable during the test.




Pain-free, pre-test planning. Brought to you by SentrySuite™.

In CPET testing, pre-test set-up can be detailed and time-consuming. SentrySuite™ provides easy to use tools for pre-exercise decision making and protocol modification. Plus, all pre-test questionnaire information can be collected via iPad and automatically uploaded to SentrySuite. Everything is right where you need it, when you need it most.


CPET Start-Up Menu: All Pre-Exercise Decisions on a Single Screen.

- Color-coded hardware connection check
- Provides suggested target load and automated protocol selection, based on measured PFT values and max predicted values
- Choose preferred test layouts, mask size, and breath averaging from start-up screen



Profile Editor Tool: Build powerful automated Protocol and Submeasurement Programs


- Easily create individual ramp, step and weight dependent protocols
- Add measurements including automated BP, RPE, exercise flow volume loops, lactates, and blood gas
- Coordinated graphical and tabular representation of events
- Multi-stepdown recovery capability




And, intuitive step-wise evaluation for occasional and frequent users.

Using step-by-step guidance, SentrySuite™ makes post-test evaluation simple and systematic. Now, evaluation and interpretation can be standardized, reducing time to result. And, workflows can be configured for individual users in relation to desired tasks and sequences. For experts, SentrySuite also provides a easy post-test way to enter offline blood gases for automatic P(A-a) O₂ and VD/VT calculation.

Workflow Steps



Accept SentrySuite's automatic selection of baseline and peak data, or manually over-ride with simple click and drag.



Side-by-side graphics with plausibility checks makes viewing ventilatory thresholds accurate and easy. Adjust the range of interest in one graph; all subsequent graphs and tabular values adjust automatically.

Focused and integrated.

Complete data review and reporting is both intuitive and automated.

- 1 Choose breath or time averaging
- 2 Quickly view, print or store reports
- 3 Quickly search for stored markers like e.g. lactate or blood gases
- 4 Quick patient-data access
- 5 Tabs to quickly switch to view different graphics
- 6 Tabular data with adjustable filtering/averaging
- 7 Start edit mode for thresholds, slopes, ranges or exercise flow volume loops (EFVL)
- 8 View/Hide recovery data from graphical displays
- 9 Comments/interpretation tool with user-definable templates and automated CPET interpretation included
- 10 Color-coded classification bar based on 'VO₂ Max predicted'
- 11 Post-test Templates and Macros

CPET BxB SENTRYSUITE V3.20.3
 Patient: Vyaire CPET
 ID: Max_bike Visit: 08/04/2022
 Status: CPET BxB Result mode

Results Table

Time min	Load W	HR BPM	O2pulse mL	Psys mmHg	Pdia mmHg	VE L/min	VTex L	VO2/kg (mL/min/kg)	SpO2 %	VO2 L/min	VCO2 L/min	RER	BR FEV% %	PETCO2 mmHg
00:21	0	61	4.3			10	0.586	4.1	97	0.3	0.2	0.76	92	29.92
00:25	0	61	4.3			10	0.573	4.1	97	0.3	0.2	0.76	92	29.64
00:29	0	61	4.4			11	0.581	4.1	97	0.3	0.2	0.77	92	29.78
00:32	0	61	4.5			11	0.587	4.2	97	0.3	0.2	0.77	92	29.64
00:34	0	61	4.6			11	0.599	4.3	97	0.3	0.2	0.77	92	29.91
00:37	0	62	4.5			11	0.597	4.3	97	0.3	0.2	0.78	92	29.72
00:40	0	62	4.5			11	0.596	4.3	97	0.3	0.2	0.78	92	29.63
00:44	0	62	4.5			11	0.599	4.3	96	0.3	0.2	0.78	92	29.66
00:47	0	62	4.4			11	0.602	4.2	96	0.3	0.2	0.78	92	29.76
00:51	0	62	4.3			10	0.612	4.1	96	0.3	0.2	0.78	92	29.88
00:55	0	62	4.4			10	0.627	4.2	96	0.3	0.2	0.78	92	30.14
00:58	0	62	4.5			10	0.622	4.2	97	0.3	0.2	0.78	92	30.19

Graphs:
 - **VO2 / Load / HR / VCO2:** Shows oxygen uptake, heart rate, and carbon dioxide output over time. Markers for VT1 and VT2 are indicated.
 - **RER / SpO2 / BR FEV%:** Shows respiratory exchange ratio, oxygen saturation, and breathing reserve over time. Markers for VT1 and VT2 are indicated.

PRE EXERCISE SPIROMETRY

	PRE	PRED	%PRED
FVC	L 4.37	4.72	92.6
FEV1	L 3.60	3.77	95.5
FEV1/FVC	% 82	80	102.8
MVV	L/min 130		

CARDIOVASCULAR RESPONSES

Summary	Rest (R) 8-Breaths	AT (VT1) [Manual] 8-Breaths	MEASURED 8-Breaths	PEAK PRED 8-Breaths	%PREDICTED 8-Breaths
VO2/kg [(mL/min)/kg]	5.0	19.0	43.0	117	117
VO2 [L/min]	0.3	1.2	2.8	117	117
VCO2 [L/min]	0.3	1.1	3.7		
Load [W]	0	100	287	159	159
HR [BPM]	64	94			90
O2pulse [mL]	5.0	13.1			131
HRR [BPM]	111	81			
Psys [mmHg]		142			


Interpretation of visit:
 CPET Eschenbacher, Mannina (1990)
 PULMONARY LIMITATION TO EXERCISE
 Mild ventilatory mechanical limitation
 CARDIAC OR CIRCULATORY LIMITATION TO EXERCISE
 If ischemic symptoms (chest pain, ST-segment changes, etc) are present:
 Ischemic heart disease
 If ischemic symptoms (chest pain, ST-segment changes, etc) are absent:
 No obvious cardiac or circulatory limitation
 This is a computer interpretation, review by a physician is required.

Classification based on VO2% Predicted:
 10 Severe Moderate Mild Normal Excellent

Physician Summary:
 A cardiopulmonary stress test was performed on CPET Vyaire using Vyntus CPX metabolic system. The patient exercised on a cycle ergometer using a ramp protocol. Heart rate was measured using the Vyntus ECG system. The results of the test are as follows:
 FUNCTIONAL CAPACITY:
 The patient exercised for 11:37 minutes. The peak work rate achieved was 287 W. The maximum oxygen uptake (VO2 Peak) was 2.8 L (mL/min)/kg, 117.4% of predicted.
 The ventilatory threshold (VT1) occurred at 1.2 L/min, 51.9% of predicted max VO2 (normal >40%).
 The ventilatory threshold (VT2) occurred at 2.8 L/min, 117.0% of predicted max VO2 (normal > 60%).
 RESPIRATORY/VENTILATORY:
 The minute ventilation (VE) at max was 108 L/min. The breathing reserve (BR = (FEV1*xx - VE_max)/FEV1*xx) was 20% (normal 20-40%).
 The VD/VT (estimated from PetCO2) at rest was 20% (normal 30%) and at maximum exercise 5% (normal less than 20%).

CARDIOPULMONARY EXERCISE DATA COMPARISON- CURRENT TO PREVIOUS

Identification: TEST 2
 Date of Birth: 5/18/1977
 Age: 41 Years
 Last Name: OD
 Height: 69.0 Inch
 Weight: 180.0 lbs
 First Name: B
 Gender: male
 BMI: 27 kg/m²



	NOW	PREV	% DIFF	PREV -1	PREV -2	PREV -3	PREV -4	PREV -5
Measurement date	10/25/18	10/25/18		03/29/16	12/16/15			
Age	41 Years	38 Years		38 Years	38 Years			

SLOPES		NOW	PREV	% DIFF	PREV -1	PREV -2
VECO2s		23.08	31.16	-25.9	31.65	27.48
VO2Ws	(mL/min)/W	5.63	10.71	-47.4	10.02	10.40

SPIOOMETRY		NOW	PREV	% DIFF	PREV -1	PREV -2
FVC	L	4.81	4.81	0.0	5.52	5.46
FEV1	L	3.85	3.85	0.0	4.24	4.23
MVV	L/min				203	189

PEAK CARDIOVASCULAR RESPONSES		NOW	PREV	% DIFF	PREV -1	PREV -2
V'O2	L/min	3.4	3.2	6.3	3.1	3.5
V'O2/kg	(mL/min)/kg	41.8	34.5	21.1	33.0	37.7
V'CO2	L/min	3.9	4.0	-2.3	3.9	4.0
Load	W	415	284	46.1	292	312
HR	BPM	171	162	5.6	171	181
O2pulse	mL	20.0	19.8	0.7	17.9	19.3
Psys	mmHg	150	187	-19.8	195	140
Pdia	mmHg	80	92	-13.1	87	86

PEAK VENTILATORY RESPONSES		NOW	PREV	% DIFF	PREV -1	PREV -2
V'E	L/min	94	123	-23.6	124	132
VTex	L	3.286	3.924	-16.3	3.677	3.709
RR	1/min	28.7	31.5	-8.8	33.7	35.5
BR FEV%	%	39	8	359.4	16	11

PEAK GAS EXCHANGE RESPONSES		NOW	PREV	% DIFF	PREV -1	PREV -2
RER		1.14	1.24	-8.0	1.26	1.14
SpO2	%	93	98	-5.1	100	100
EqCO2		23.2	30.1	-22.8	31.0	31.8
EqO2		26.4	37.3	-29.0	39.1	36.3
PETCO2	mmHg	47.13	36.17	30.3	35.89	35.17
PETO2	mmHg	103.31	117.10	-11.8	117.50	117.59

Trending patient data over time is key.

SentrySuite™'s ability to compare a patient's metabolic data longitudinally, in both graphical and tabular forms, has taken on a whole new meaning in the post COVID-19 era. SentrySuite's unique ability is more powerful than ever.

Report Generation, our unique report designer, is both simple to use, and extremely comprehensive.

With SentrySuite™, eye-catching results available on every page.




Of course SentrySuite provides all standard reports, including the 9-panel graph. But where it really shines is its capacity to redefine how you visually present your data in more powerful and meaningful ways. Three great examples include our Exercise Flow Volume Loops, CPET comparison graph, and Dynamic Predicted Ranges.

Exercise Flow Volume Loops

Quickly see dynamic airway hyperinflation and flow limitation by trending Flow/Volume and EILV/EELV throughout exercise.

Dynamic Predicted Ranges

Continually monitor VE/VCO₂ data during test.



Graph 1 shows a normal response to exercise and graph 2 shows an elevated response to exercise

Wide-ranging data viewing including easy selection of steady state data.

The ability to easily and quickly trend data is equally important when performing ongoing, routine QA protocols to assess the health status of your device. Here we are showing a Bio-QC bike 3-stage steady state protocol that is standard to SentrySuite™ profiles.



Resting Energy Expenditure

Simple mask methodology; no new hardware

Easily view when patient reaches steady-state conditions

Select up to four areas of steady-state showing data averages with coefficient of variation (CV)

Steady State Qualitative Analysis

Perform monthly as part of routine QA plan

3 simple bike steps under AT or use any published treadmill QA protocol

Average last 2 minutes of each step and compare to published norms

Expand your capability by combining Vyntus™ CPX with other devices.



GE CardioSoft® ECG



Tango® blood pressure monitor



Polar® Bluetooth® interface



Ergoselect 600P recumbent bike



Treadmill



VIAsprint® 150 / 200P cycle ergometer with / without BP

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA

REG ONLY For U.S. distribution only

Trademarks are the property of their respective owners.
© 2022 Vyaire. Vyaire, the Vyaire logo and all other trademarks or registered trademarks are property of Vyaire Medical, Inc., or one of its affiliates. VYR-US-1800020 1.0

