



Respiratory ventilation products

Portable ventilators

Table of contents

LTV™ Series	3	LTV Series miscellaneous parts & accessories	20-21
LTV Series ventilators	4	LTV Series ventilator miscellaneous parts & accessories	22
LTV Series accessories.....	5	LTV Series accessories, manuals & educational materials.....	23
LTV Series circuits.....	6-10	ReVel™ ventilator	24
LTV Series circuit accessories.....	11	ReVel ventilator accessories.....	25-29
LTV Series power accessories and batteries	12-15	Terms and Conditions	30
LTV Series ventilator stands, packs and mounting accessories.....	16	Index	31
LTV Series stands, packs and mounting accessories.....	17-19		

LTV™ Series

LTV Series ventilators	4
LTV Series accessories.....	5
LTV Series circuits.....	6-10
LTV Series circuit accessories.....	11
LTV Series power accessories and batteries	12-15
LTV Series ventilator stands, packs and mounting accessories	16
LTV Series stands, packs and mounting accessories.....	17-19
LTV Series miscellaneous parts & accessories.....	20-21
LTV Series ventilator miscellaneous parts & accessories	22
LTV Series accessories, manuals & educational materials.....	23





21102-001

LTV Series 1100 ventilator system

LTV Series 1100 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter and operator manual



18984-001

LTV Series 1150 system

LTV Series 1150 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter, operator manual and quick reference guide



18888-001

LTV Series 1200 system

LTV Series 1200 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter, low flow oxygen adapter, operator manual and quick reference guide



18888-201

LTV Series 1200 MR Conditional ventilator system

LTV 1200 MR ventilator with protective boots, AC adapter with power cord, dust cover, circuit adapter, low flow oxygen adapter, LTV 1200 MR Conditional floor stand, case of 5 MR-safe 15 ft circuits and operator manual



99503


LTV 1200 MR Conditional System with SprintPack™

LTV 1200 MR Conditional ventilator, protective boots, AC adapter with power cord, SprintPack™ bracket, SprintPack™ battery power system, dust cover, circuit adapter, low flow oxygen adapter, LTV 1200 MR Conditional floor stand, case of 5 MR-safe 15 ft circuits and operator manual

LTM™ II Graphics Monitor and accessories



17650-001
LTM II Graphics Monitor package
 LTM II Graphics Monitor, LTM mounting bracket, LTM II data cable, power splitter cable, compact flash card and LTM II operator manual



17916-001
Bracket, LTM Graphics Monitor mounting



11089
Cable, LTM data, 21 in



11089-002
Cable, LTM data, 6 ft



18093-001
Cable, power splitter, LTM II




18097-001
Stylus pen, LTM



18098-001
Card, compact flash

LTV Series 1100, 1150 and 1200 ventilator circuits without PEEP


Single patient use (SPU) circuits, 6 ft



29657-001

Patient circuit without PEEP, 22 mm SPU (10/pkg)


1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb



29658-001

Patient circuit without PEEP, 15 mm SPU (10/pkg)


1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb



29663-001

Patient circuit without PEEP with 1 water trap, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb with water trap, temp probe port, wye and 61 cm x 22 mm exp limb




29662-001

Patient circuit without PEEP with 1 water trap, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb with water trap, temp probe port, wye and 61 cm x 15 mm exp limb


Reusable circuits, 6 ft



17997-001

Patient circuit without PEEP, 22 mm reusable (ea) (replaces 10815)

1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb





18941-001

Patient circuit without PEEP, 15 mm reusable (ea)

1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb

LTV Series 1100, 1150 and 1200 ventilator circuits without PEEP (continued)


SPU heated wire circuits for Fisher & Paykel humidifiers, 5.5 ft



15090-102

Patient circuit without PEEP, 22 mm SPU single heated wire, Fisher & Paykel (10/pkg)



1.7 m x 22 mm heated insp limb, temp probe port, wye and 61 cm x 22 mm exp limb



15090-104

Patient circuit without PEEP, 15 mm SPU single heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye and 61 cm x 15 mm exp limb



15091-102

Patient circuit without PEEP with elbow, 22 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 22 mm heated insp limb, temp probe port, wye, elbow and 61 cm x 22 mm heated exp limb




15091-104

Patient circuit without PEEP with elbow, 15 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye, elbow and 61 cm x 15 mm heated exp limb

SPU MR safe circuit, 15 ft




19189-001

Patient circuit without PEEP, 22 mm SPU, 15 ft (5/pkg)

LTV Series 900, 950 and 1000 ventilator circuits with PEEP



SPU circuits, 6 ft



29695-001

Patient circuit with PEEP, 22 mm SPU (10/pkg)



1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb



29710-001

Patient circuit with PEEP, 15 mm SPU (10/pkg)



1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb



29697-001

Patient circuit with PEEP with 1 water trap, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb with water trap, temp probe port, wye and 61 cm x 22 mm exp limb





29696-001

Patient circuit with PEEP with 1 water trap, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb with water trap, temp probe port, wye and 61 cm x 15 mm exp limb

LTV Series 900, 950 and 1000 ventilator circuits with PEEP (continued)


SPU heated wire circuits for Fisher & Paykel humidifiers, 5.5 ft



15090-101

Patient circuit with PEEP, 22 mm SPU single heated wire, Fisher & Paykel (10/pkg)



1.7 m x 22 mm heated insp limb, temp probe port, wye and 61 cm x 22 mm exp limb



15090-103

Patient circuit with PEEP, 15 mm SPU single heated wire, Fisher & Paykel (10/pkg)


1.7 m x 15 mm heated insp limb, temp probe port, wye and 61 cm x 15 mm exp limb



15091-101

Patient circuit with PEEP, 22 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 22 mm heated insp limb, temp probe port, wye and 61 cm x 22 mm heated exp limb




15091-103

Patient circuit with PEEP, 15 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye and 61 cm x 15 mm heated exp limb

LTV Series 800 ventilator circuits with PEEP and without flow sensor (for pressure triggering)


SPU circuits, 6 ft



29715-001

Patient circuit with PEEP, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb, wye and 61 cm x 22 mm exp limb



29716-001

Patient circuit with PEEP, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb, wye and 61 cm x 15 mm exp limb



10696
Adapter, circuit



17921-001
Exhalation valve with PEEP, reusable
(replaces 10824)



18172-001
Exhalation valve without PEEP, reusable
(replaces 10825)



17486-101
Diaphragm, exhalation valve
(for 17921-001 and 18172-001)

SprintPack™ Lithium-Ion Power System



19222-201
 Power system, SprintPack lithium-ion with cable protector
 SprintPack power manager, cable protector, 2 SprintPack batteries and operator manual



19333-201
 Power manager, SprintPack with cable protector (Does not include SprintPack batteries)



19444-001
 Battery, SprintPack



19096-001
 Bracket, SprintPack mounting



19102-001
 Transport pack for LTV and SprintPack
 5 D-rings, detachable shoulder strap, hand strap and bedrail mount



11544
 Cord, automotive power, 12 ft



19871-001
 Mounting bracket, MR Conditional SprintPack



24555-001
 SprintPack Cable Support Bracket

External sealed lead-acid batteries



10800
 Battery kit, 9-hour external
 12 V/33 amp-hour battery, battery box with strap, charger, battery cable and operator manual



10787
 Battery, 9-hour external (12 V/33 amp-hour)



10790
 Battery box with strap, external



10900
 Battery kit, 3-hour external
 12 V/12 amp-hour battery, padded battery bag with cable, charger, battery cable and operator manual



11475
 Battery, 3-hour external (12 V/12 amp-hour)



10901
 Cable, 3-hour external battery



11330
 Bag with cable, padded battery
 (for use with 3-hour battery)



10802
 Cable, DC external battery
 (for use with 3-hour and 9-hour battery)



10801
 Charger, external battery
 (for use with 3-hour and 9-hour battery)

Universal power supply (UPS)



17894-001
Battery, universal power supply

Miscellaneous power accessories



11448
AC power adapter with power cord, LTV



17527-001
AC power adapter without power cord, LTV



10854
Fuse, 32 V, 10 amp battery



11456
Electrical strip with 4 m power cord, 4-outlet U.S. hospital grade



11544
Cord, automotive power, 12 ft



10536
Power cord, 110 VAC hospital grade, 10 ft

LTV Series stands



10611
Stand with oxygen cylinder ring,
LTV Series floor



10695
Stand, LTV Series table

LTV Series mounts, packs and carriers




10700
Mount, bedrail



11581-G
Backpack for LTV Series ventilators with boots *(gray)*



19102-001
Transport pack for LTV Series ventilators and SprintPack
5 D-rings, detachable shoulder strap, hand strap and bedrail mount



19184-001
Freedom carrier, LTV Series ventilators and SprintPack

LTV Series mounting accessories




10610
Bracket, female dove-tail



10979
Bracket, humidifier pole mount C-clamp



11454
Bracket, large battery mounting



17916-001
Bracket, LTM graphics monitor mounting



19096-001
Bracket, SprintPack mounting



10658
Hand strap, LTV Series ventilators
(gray)



19031-001
Hand strap, LTV Series ventilators
(purple)



11501
Support arm with circuit hanger, patient circuit

LTV Series mounting accessories (continued)



11503
Hanger, clip, circuit support arm



11450
Basket, utility



11451
Crossbar, short



11452
Crossbar, long



11453
Pole, humidifer



11455
Oxygen cylinder ring set
(for use with 10611)



17341-001
Strap, transport pack shoulder

LTV Series mounting accessories sold but not shown

Cat. no.	Description
CAS-1152	Caster for floor stand (for use with 10611)
CAS-1168	Caster for floor stand, locking (for use with 10611)

LTV Series mounting accessories (continued)



19870-001
L-bracket with screws for LTV Series MR Conditional floor stand



19871-001
Mounting bracket, LTV Series MR Conditional SprintPack



19873-001
Caster for LTV Series MR Conditional floor stand



19874-001
Caster for LTV Series MR Conditional floor stand, locking



19875-001
Oxygen cylinder ring set, LTV Series MR Conditional floor stand



34910
Cable, remote alarm, 50 ft
(for use with alarm 11553)



19103-001
Cable, remote alarm, BNC - phone, 20 ft



14647-001
Cable, nurse call splitter



10779
Cable, nurse call, normally closed



10780
Cable, nurse call, normally open



17986-001
Adapter, electrical, Fisher & Paykel 900 MR
900 dual heated limb




17987-001
Adapter, electrical, Fisher & Paykel 900 MR
901 single heated limb



10701
Adapter, low pressure oxygen



10702
Adapter, NIST oxygen




11550
Boot kit
(gray protective)



19175-001
Boot kit
(purple protective)



15016-001
Cap, blender



14446
Cap, tethered port
(for use with LTV Series 800, 900, 950, 1100 and 1150 ventilators)



18098-001
Card, compact flash




10109
Cover, front panel dust




11456
Electrical strip with 4 m power cord,
4-outlet U.S. hospital grade



10678
Fan grill assembly
Fan grill, fan grill enclosure, fan filter



10258x5
Filter, LTV Series air inlet (5/pkg)



10789x5
Filter, LTV Series fan (5/pkg)




10699
Hose, 15 ft oxygen
(green)



11466
Test lung, 1L



10934x50
Tie, circuit twist-loc, dual (50/pkg)



10936x50
Tie, circuit twist-loc, large (50/pkg)

LTV Series ventilator shipping carton



Cat. no.	Description
19818-001	Carton, shipping (replaces 17496-001)
17490-001	Sheet, divider (secures AC adapter)
17489-001	Sheet, divider (secures operator manual)
17487-001	Cap set, LTV end

Order all four parts

Cat. no.	Description
19802-001	Operator's Manual, LTV 1200/1150 English
19822-001	Quick Reference Guide, LTV 1200/1150 English
18999-001	In-Service DVD, LTV 1200 English
19990-001	Caregiver Training DVD, LTV 1150 English
19960-001	Emergency Set-Up DVD, LTV 1200 English
21106-001	Operator's Manual, LTV 1100 English
19795-001	Operator's Manual, LTV 1000/950/900 English
19819-001	Quick Reference Guide, LTV 1000/950/900 English
18581-001	In-Service CD/DVD, LTV 1000/950/900 English
18580-001	Caregiver Training DVD, LTV 1000/950/900 English
11600	Simulator CD, LTV 1000 English
17622-001	Recommended Settings Card, LTV 1000/950/900 English (25/pkg)
19797-001	Operator's Manual, LTV 800 English
19821-001	Quick Reference Guide, LTV 800 English
11524	Quick Reference Guide, HME/Hospital English
17754-001	Operator's Manual, LTM II Graphics Monitor English
11009	LTM Graphics Monitor CD, English
21188-001	Operator's Manual, LTV 1200/1150 Spanish
21196-001	Quick Reference Guide, LTV 1200/1150 Spanish
19799-001	Operator's Manual, LTV 1000/950/900 Spanish
19820-001	Quick Reference Guide, LTV 1000/950/900 Spanish

ReVel™ ventilator

ReVel ventilator accessories.....25-30



Docking stations, stands and accessories



21380-001
Protective collar, docking station



13914-001
Basket, O₂ cylinder (for use with floor stand 13912-001)



19574-001
Stand, table



13886-001
Wall mount, humidifier pole



13912-001
Stand, floor, with patient circuit arm
Floor stand, patient circuit support arm, 3/16" hex wrench and instructions



12140-001
Arm, patient circuit support

Docking stations, stands and accessories *(continued)*



19937-001
Mounting plate, docking station
(for use with GCX wall mounts)



12695-002
ReVel carry case



24737-001
PTV Transport Bracket for bed rail and pole mounting. Includes PalmTop mounting accessory and clamp.
Packaged preassembled.



24736-001
PTV Transport Bracket for L-Track mounting. Includes PalmTop mounting accessory and semi-permanent hardware.
Assembly required.



18098-001
Card, compact flash

Power accessories and batteries



12103-001
Power adapter with power cord, 110 VAC



13955-001
Charger, desktop battery with 110 VAC power cord *(batteries not included)*



10536
Power cord, 110 VAC, hospital grade, 10'

Power accessories and batteries *(continued)*



13908-001
Battery, main lithium-ion



12102-001
Adapter, automobile power point, 12'

Miscellaneous parts and accessories



10699
Hose, O₂, 15'




13036-001
Sensor, O₂



13037-001
T-Piece, reusable O₂ sensor



13038-001
T-Piece, SPU O₂ sensor (10/pkg)



13897-001
Cable, O₂ sensor



18602-001
Cable, SpO₂

Miscellaneous parts and accessories (continued)



13063-001
Adapter and cable, VOXP data



19826-001
Filter, air inlet



19833-001
Filter, fan (10/pkg)




11466
Test lung, 1L

ReVel ventilator accessories sold but not shown

Cat. no.	Description	Qty.
19463-001	SpO ₂ Flexi-Form II, SPU adult	10/pkg
19464-001	SpO ₂ Flexi-Form II, SPU pediatric	10/pkg
19465-001	SpO ₂ Flexi-Form II, SPU infant	10/pkg
21311-001	SpO ₂ Flexi-Form III, SPU adult	24/pkg

ReVel circuits


SPU circuits, 6'



14858-001

Patient circuit without PEEP with elbow, 22mm SPU (10/pkg)

with 1.8 m x 22 mm ID inspiratory limb, wye, elbow and 61 cm x 22 mm expiratory limb




14859-001

Patient circuit without PEEP with elbow, 15mm SPU (10/pkg)

with 1.8 m x 15 mm ID inspiratory limb, wye, elbow and 61 cm x 15 mm expiratory limb


SPU heated wire circuits for Fisher & Paykel (F&P) humidifiers, 5.5'



15091-102

Patient circuit without PEEP with elbow, 22 mm SPU dual heated wire, F&P (10/pkg)

with 1.7 m x 22 mm ID heated inspiratory limb, temp probe port, wye, elbow and 61 cm x 22 mm heated expiratory limb



15091-104

Patient circuit without PEEP with elbow, 15 mm SPU dual heated wire, F&P (10/pkg)

with 1.7 m x 15 mm ID heated inspiratory limb, temp probe port, wye, elbow and 61 cm x 15 mm heated expiratory limb

Circuit accessories

Cat. no.	Description
18225-001	Patient wye and sense lines, reusable

Minimum order of \$50.00 required

Terms and Conditions of Sale:

Pricing and warranties are subject to change without notice.

Merchandise Return Policy:

1. Please contact our Customer Care Team for a Return Material Authorization (RMA) number prior to returning merchandise
2. Customers may return unopened and unused product within 30 days
3. All products are subject to a 20% restocking fee
4. Please be prepared with the following information:
Item number, Serial or Lot Number, with Purchase Order or Sales Order Number

To place an order or for more information, please contact Customer Service at:

Phone:

+1 800-754-1914

Fax:

+1 763-398-8359

Email:

rcorders@vyaire.com

Numerical index

LTV Series ventilators						ReVel ventilators	
10109	21	11544	14	19102-001	16	10536	26
10258x5	21	11550	21	19103-001	20	10699	27
10536	14	11581-G	16	19175-001	21	11466	28
10610	17	11600	23	19184-001	16	12102-001	27
10611	15	14446	21	19189-001	7	12103-001	26
10658	17	14647-001	20	19222-201	12	12140-001	25
10678	21	15016-001	21	19333-201	12	12695-002	26
10695	15	15090-101	9	19444-001	12	13036-001	27
10696	11	15090-102	7	19795-001	23	13037-001	27
10699	22	15090-103	9	19797-001	23	13038-001	27
10700	16	15090-104	7	19799-001	23	13063-001	28
10701	20	15091-101	9	19802-001	23	13886-001	25
10702	20	15091-102	7	19818-001	23	13897-001	27
10779	20	15091-103	9	19819-001	23	13908-001	27
10780	20	15091-104	7	19820-001	23	13912-001	25
10787	13	17341-001	18	19821-001	23	13914-001	25
10789x5	22	17486-101	11	19822-001	23	13955-001	26
10790	13	17487-001	23	19870-001	19	14858-001	29
10800	13	17489-001	23	19871-001	12	14859-001	29
10801	13	17490-001	23	19871-001	19	15091-102	29
10802	13	17527-001	14	19873-001	19	15091-104	29
10854	14	17622-001	23	19874-001	19	18098-001	26
10900	13	17650-001	5	19875-001	19	18225-001	29
10901	13	17754-001	23	19960-001	23	18602-001	27
10934x50	22	17894-001	14	19990-001	23	19463-001	28
10936x50	22	17916-001	5	21102-001	4	19464-001	28
10979	17	17916-001	17	21106-001	23	19465-001	28
11009	23	17921-001	11	21188-001	23	19574-001	25
11089	5	17986-001	20	21196-001	23	19826-001	28
11089-002	5	17987-001	20	24555-001	12	19833-001	28
11330	13	17997-001	6	29657-001	6	19937-001	26
11448	14	18093-001	5	29658-001	6	21311-001	28
11450	18	18097-001	5	29662-001	6	21380-001	25
11451	18	18098-001	5	29663-001	6	24736-001	26
11452	18	18098-001	21	29695-001	8	24737-001	26
11453	18	18172-001	11	29696-001	8		
11454	17	18580-001	23	29697-001	8		
11455	18	18581-001	23	29710-001	8		
11456	14	18888-001	4	29715-001	10		
11456	21	18888-201	4	29716-001	10		
11466	22	18941-001	6	34910	20		
11475	13	18984-001	4	99503	4		
11501	17	18999-001	23	CAS-1152	18		
11503	18	19031-001	17	CAS-1168	18		
11524	23	19096-001	12				
11544	12	19096-001	17				
		19102-001	12				



GLOBAL HEADQUARTERS

Vyairē Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045, USA

[vyairē.com](https://www.vyairē.com)

 **For U.S. distribution only.**

©2020 Vyairē. Vyairē and the Vyairē Logo are trademarks of Vyairē Medical, Inc. VYR-US-2000005