



# Disposables and Accessories Catalog

Pulmonary Function Testing  
Cardiopulmonary Exercise Testing

# About this catalog

**We are pleased that you have chosen to use Vyair  
Medical products.**

The Vyair product portfolio for Respiratory Diagnostics ranges from small mobile devices to large pulmonary function laboratory devices. As a result there is a variety of consumables and accessories associated with each specific product. This catalog is designed to ensure that you find the correct consumables and accessories for the specific piece of Vyair equipment you are using both quickly and effectively.

**For ease of use the catalog has been divided into three hierarchical levels:**

- 1. Level 1, functional areas such as Pulmonary Function Testing, Cardiopulmonary Exercise Testing, etc.**
- 2. Level 2, subdivides the functional areas into specific device names. Look for the name of your device. If the device name is not very clear, the device overview on the following pages will give you assistance.**
- 3. Level 3, is a further subdivision of the devices into their measuring applications or product groups e.g. Spirometry.**

An index directory at the end of the catalog can also be used to search for all consumables and accessories listed in the catalog.

Vyair consumables and accessories are medical products and are designed, selected and validated specially for Vyair devices. Please use only the original consumables and accessories from Vyair to ensure accuracy of results and patient safety.

We are happy to hear from you.  
Your Vyair Team



**Our phone order service:**

**USA:**  
1.800.231.2466 Toll-Free  
1.714.283.2228 Phone



**Mail:**

**USA:**  
RCorders@vyair.com



# Table of Contents

## Pulmonary Function Testing (PFT) 8

### Vyntus™ BODY | Vyntus™ ONE | Vyntus™ SPIRO

MicroGard™ II Bacterial/Viral Filter .....	10
Spirometry, Body Plethysmography, Diffusion, N <sub>2</sub> Washout .....	11
Gases and Cylinders for Diffusion and N <sub>2</sub> Washout .....	19
Volume Calibration.....	20

## Cardiopulmonary Exercise Testing (CPET), ECG and 6-Minute-Walk-Test 22

### Vyntus™ CPX | Vyntus™ ONE

Masks, Mouthpieces and Nose Clips .....	24
SpO <sub>2</sub> Sensors .....	28
Calibration Gases.....	31
Standard Breath by Breath.....	32
High/Low FIO <sub>2</sub> .....	34

### Vyntus™ ECG

All Accessories .....	36
-----------------------	----

### CardioSoft EKG CAM-USB

All Accessories .....	37
-----------------------	----

### Blood Pressure Accessories for VIA sprint | Ergoselect Series | Suntech Tango

All Accessories .....	38
-----------------------	----

### Ergometer Accessories for VIA sprint 150P

All Accessories .....	40
-----------------------	----

### Vyntus™ WALK

All Accessories .....	42
-----------------------	----



# Device Overview

## Pulmonary Function Testing (PFT)



Vyntus™ SPIRO



Vyntus™ ONE



Vyntus™ BODY

## Cardiopulmonary Exercise Testing (CPET), ECG and 6-Minute-Walk-Test



Vyntus™ CPX



Vyntus™ ONE



Vyntus™ ECG



GE CardioSoft EKG  
CAM-USB / CAM-14



VIAsprint 150P with BP



Ergoselect 600 and  
600P Recumbent with  
BP



Tango M2



Vyntus™ WALK



# Pulmonary Function Testing (PFT)

Vyntus™ BODY, Vyntus™ ONE, Vyntus™ SPIRO




The following parts are suitable for these products:

Vyntus™ BODY | Vyntus™ ONE | Vyntus™ SPIRO

**MicroGard™ II Bacterial/Viral Filter**

<b>MicroGard™ IIB single use bacteria filter with integrated mouthpiece</b>	For Device	
	All PFT devices	

Package	50 pieces
Package	V-892381

<b>Filter Kit MicroGard™ Type IIB single use</b>	For Device	
	All PFT devices	

Package	80 pieces	Contains one filter and one nose clip in individually wrapped package. Delivered in box of 80.
Package	V-892391	

<b>MicroGard™ IIC single use bacteria filter</b>	For Device	
	All PFT devices	

Package	50 pieces	Can be used in combination with mouthpieces. Connection to mouthpiece outer diameter 30 mm.
Package	V-892384	

<b>Filter Kit MicroGard™ Type IIC single use</b>	For Device	
	All PFT devices	

Package	80 pieces	Contains one filter, one mouthpiece and one nose clip in individually wrapped package. Delivered in box of 80.
Package	V-892392	


**Spirometry, Body Plethysmography, Diffusion, N<sub>2</sub> Washout**

<b>Nose clip, single use</b>	For Device	
	All PFT devices	

Package	10 pieces
Part Number	42101

<b>Nose clip, single use</b>	For Device	
	All PFT devices	


Package	100 pieces
Part Number	761889


<b>Mouthpiece plastic, single use</b>	For Device	
	All PFT devices	

Package	25 pieces
Part Number	001-892103


<b>Mouthpiece silicone adult, latex free, reusable</b>	For Device	
	All PFT devices	


Package	5 pieces
Part Number	V-892235

<b>Mouthpiece silicone children, latex free, reusable</b>		For Device	
		All PFT devices	
Package	5 pieces		
Package	V-892230		


<b>Mouthpiece ID30, single use</b>		For Device	
		All PFT devices	
Package	180 pieces		
Package	V-892894		


<b>Pneumotach plastic PT 36 complete, not heatable</b>		For Device	
		Vyntus SPIRO	
Package	1 piece		
Package	001-709554		


<b>Vyntus™ SPIRO pneumotach/ elbow kit</b>		For Device	
		Vyntus SPIRO	
Package	1 kit	Contains four pneumotachs (001-709554) and five elbows (V-892650).	
Package	33120-001		


<b>Spare screen PT 36</b>		For Device	
		Vyntus SPIRO	
Package	1 piece		
Part Number	001-700096		


<b>O-ring for plastic P-tube PT 36, D 63 x 2 mm</b>		For Device	
		Vyntus SPIRO	
Package	1 piece		
Part Number	001-851243		


<b>Elbow 45°, ID 30 mm</b>		For Device	
		Vyntus SPIRO	
Package	5 pieces		
Part Number	V-892650		


<b>Elbow 45°, ID 30 mm, ULTEM, especially resistant against disinfection</b>		For Device	
		Vyntus SPIRO	
Package	1 piece		
Part Number	001-852360		


<b>Sample line Vyntus™ Body 48 inch</b>		For Device	
		Vyntus BODY	
Package	1 piece		
Package	V-852556		


<b>Sample line Vyntus™ ONE 240 cm</b>		For Device	
		Vyntus ONE	
Package	1 piece		
Package	V-707722		


<b>USS Module with adapter and one protective cover</b>		For Device	
		Vyntus BODY Vyntus ONE	
Package	1 set		
Package	V-707377		

<b>Vyntus™ ONE patient CT kit</b>		For Device	
		Vyntus ONE	
Package	1 kit	Consists of 1 x V-707377 USS Module, 1 x V-707721 FPV block eDemand, 1 x V-707744 FPV shutter mechanism C-shape, 1 x V-707390 FPV silicone valves	
Package	32924-001		


<b>Vyntus™ BODY patient CT kit</b>		For Device	
		Vyntus BODY	
Package	1 kit	Consists of 1 x V-707377 USS Module, 1 x V-707389 FPV block, 1 x V-707388 FPV shutter mechanism, 1 x V-707390 FPV silicone valves	
Part Number	32925-001		


<b>Cable port protective cover</b>		For Device	
		Vyntus BODY Vyntus ONE	
Package	5 pieces		
Part Number	V-707384		

<b>USS MicroGard™ adapter</b>		For Device	
		Vyntus BODY Vyntus ONE	
Package	1 piece		
Part Number	V-707380		


<b>USS MicroGard™ adapter for gas sampling</b>		For Device	
		Vyntus BODY Vyntus ONE	
Package	1 piece		
Part Number	V-707381		





<b>FPV block</b>	For Device		
	Vyntus BODY		
Package	1 piece	Inclusive 2 x FPV silicone valve and 1 x V-707397 FPV X-Ring Kit	
Package	V-707389		


<b>FPV block eDemand</b>	For Device		
	Vyntus ONE		
Package	1 piece	Inclusive 2 x FPV silicone valve	
Package	V-707721		


<b>FPV shutter mechanism</b>	For Device		
	Vyntus BODY		
Package	1 piece		
Package	V-707388		


<b>FPV shutter mechanism C-shape</b>	For Device		
	Vyntus ONE		
Package	1 piece		
Package	V-707744		

<b>FPV silicone valves</b>	For Device		
	Vyntus BODY Vyntus ONE		
Package	5 pieces		
Part Number	V-707390		


<b>FPV housing clamp</b>	For Device		
	Vyntus BODY		
Package	1 piece		
Part Number	V-707392		


<b>FPV housing clamp eDemand</b>	For Device		
	Vyntus ONE		
Package	1 piece		
Part Number	V-861935		

<b>FPV X-Ring Kit consists of 3 x X-ring and 1 x small rod</b>	For Device		
	Vyntus BODY		
Package	1 set		
Part Number	V-707397		

<b>O<sub>2</sub>-Fuel Cell</b>		For Device
		Vyntus ONE
Package	1 piece	
Package	923015	


**Gases and Cylinders for Diffusion and N<sub>2</sub> Washout**


<b>Gas DLCO 0.3% CO, 0.3% CH<sub>2</sub> 21% O<sub>2</sub>, BAL N<sub>2</sub> (ED)</b>		For Device	
		Vyntus BODY Vyntus ONE	
Package	1 filling	•	
Part Number	768093		


<b>Gas - 99.00% OXYGEN (USP)</b>		For Device	
		Vyntus ONE	
Package	1 filling	•	
Part Number	761819		


<b>Vyntus™ ONE calibration gas 16% O<sub>2</sub>, 5% CO<sub>2</sub>, 0.3% CH<sub>2</sub>, 0.3% CO, bal N<sub>2</sub> E cylinder</b>		For Device	
		Vyntus ONE	
Package	1 filling	•	
Part Number	V-892590-13		

**Volume Calibration**

<b>Manual calibration syringe 3 L incl. silicone sleeve and certificate</b>		For Device	
		All PFT devices	
Package	1 piece		
Package	720254		

<b>QM syringe 3 L with integrated ventilator</b>		For Device	
		Vyntus BODY Vyntus ONE	
Package	1 piece		
Package	V-720000		

<b>Silicone sleeve for syringe, D 28 mm, L 60 mm</b>		For Device	
		All PFT devices	
Package	1 piece		
Package	001-852740		

<b>Silicone adapter MicroGard™ IIB to calibration syringe</b>		For Device	
		All PFT devices	
Package	1 piece		
Package	V-861449		

# Cardiopulmonary Exercise Testing (CPET), ECG and 6-Minute-Walk-Test

Vyntus™ CPX, Vyntus™ ONE,  
Vyntus™ ECG, CardioSoft ECG,  
Ergometer accessories, Blood pressure cuffs,  
Vyntus™ WALK











The following parts are suitable for these products:


Vyntus™ CPX | Vyntus™ ONE


**Masks, Mouthpieces and Nose Clips**

<b>Mask 7450 V2, reusable, size L, complete with adapter and headgear</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	V-982186		
<b>Mask 7450 V2, reusable, size M, complete with adapter and headgear</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	V-982185		
<b>Mask 7450 V2, reusable, size S, complete with adapter and headgear</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	V-982184		
<b>Mask 7450 V2, reusable, size XS, complete with adapter and headgear</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	V-982187		


<b>Mask 7450 V2, reusable, size Petite, complete with adapter and headgear</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	V-982188		


<b>Mask adapter, flexible for CPET masks 7450, 7400 and children mask 982078</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	V-982192		


<b>Headgear, size L for mask 7450 with size L</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	V-982181		


<b>Headgear, size M for mask 7450 with size M and S</b>		For Device Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	V-982180		


<b>Headgear, size S for mask 7450 with size XS and Petite</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	V-982189		

<b>Mask, Pediatric Large with head strap</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	982078		


<b>Mask, Pediatric Small with head strap</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	27934-001		

<b>Mouthpiece ID30, single use</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	180 pieces	•	
Package	V-892894		

<b>Nose clip, single use</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	10 pieces	•	
Part Number	42101		


<b>Nose clip, single use</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	100 pieces	•	
Part Number	761889		


**SpO<sub>2</sub> Sensors**


<b>Earclip SpO<sub>2</sub> sensor with 1 m cable for XPOD</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	28896-001		


<b>Extension cable SpO<sub>2</sub>, length 3 m</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Package	918321		


<b>Flex SpO<sub>2</sub> sensor, adult, length 3 m including 25 wraps</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 set	•	
Package	29836-001		

<b>Wraps for Adult Flex SpO<sub>2</sub> sensor, 8000JFW, disposable</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	25 pieces	•	
Package	22012-001		

<b>SpO<sub>2</sub> reflect sensor holder consists of 10 holder and 20 adhesive collars</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 set	•	
Part Number	28898-001		

<b>Reusable reflect sensor</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	28897-001		


<b>SpO<sub>2</sub> Finger clip sensor, Adult</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	29836-001		


<b>SpO<sub>2</sub> Finger clip sensor, Peds</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	•	
Part Number	29837-001		


<b>SpO<sub>2</sub> Kit with sensors</b>		For Device
		Vyntus CPX Vyntus ONE
Package	1 set	*Consists of 1 x V-377503 SpO <sub>2</sub> kit with ear clip, 1 x 28897-001 Reusable reflectance sensor, 1 x 28898-001 SpO <sub>2</sub> reflectance sensor holder, 1 x 29836-001 SpO <sub>2</sub> finger clip sensor adult 3m, 1 x 29837-001 SpO <sub>2</sub> finger clip sensor peds 3m
Package	29835-001	



**Calibration Gases**


<b>Calibration gas for CPET, 16% O<sub>2</sub>, 4% CO<sub>2</sub>, balance N<sub>2</sub>, 8d cylinder</b>		For Device	
		Vyntus CPX	
Package	1 filling	•	
Part Number	773867		


<b>Calibration gas for CPET, 16% O<sub>2</sub>, 4% CO<sub>2</sub>, balance N<sub>2</sub>, E cylinder</b>		For Device	
		Vyntus CPX	
Package	1 filling	•	
Part Number	673666		


<b>Vyntus™ ONE calibration gas 16% O<sub>2</sub>, 5% CO<sub>2</sub>, 0.3% CH<sub>2</sub>, 0.3% CO, bal N<sub>2</sub> E cylinder</b>		For Device	
		Vyntus ONE	
Package	1 filling	•	
Part Number	V-892590-13		





Standard Breath by Breath


<b>DVT Vyntus™ CPX/ONE, volume sensor complete</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece	Included is also one volume sensor insert 001-982094. DVT = Digital Volume Transducer.	
Package	V-707326		


<b>LED holder for DVT Vyntus™ CPX/ONE</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece		
Package	V-707325		

<b>DVT adapter for Vyntus™ CPX/ONE</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece		
Package	V-861078		


<b>TwinTube sample line 240 cm for Vyntus™ CPX</b>		For Device	
		Vyntus CPX	
Package	1 piece		
Package	V-707327		


<b>Volume sensor insert</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece		
Part Number	001-982094		


<b>Volume sensor insert, single-use</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	50 piece		
Part Number	V-707746		

<b>O<sub>2</sub>-Fuel Cell</b>		For Device	
		Vyntus CPX Vyntus ONE	
Package	1 piece		
Part Number	923015		


High/Low FIO<sub>2</sub>


<b>DVT High/Low FIO<sub>2</sub> complete for Vyntus™ CPX</b>		For Device	
		Vyntus CPX	
Package	1 piece		
Package	V-707332		


<b>LED holder High/Low FIO<sub>2</sub> for DVT Vyntus™ CPX</b>		For Device	
		Vyntus CPX	
Package	1 piece		
Package	V-707329		

<b>Y-Valve 2730 HR for Vyntus™ CPX</b>		For Device	
		Vyntus CPX	
Package	1 piece		
Package	V-861109		

<b>Adapter HR valve to breathing bag</b>		For Device	
		Vyntus CPX	
Package	1 piece		
Package	707146		

<b>Breathing bag 20 L transparent, 2 fittings</b>		For Device	
		Vyntus CPX	
Package	30 pieces		
Part Number	001-852828		

<b>Filling assembly for breathing bag</b>		For Device	
		Vyntus CPX	
Package	1 set		
Part Number	706823		

<b>Head Support for Hans Rudolph Y-Valve 2730</b>		For Device	
		Vyntus CPX	
Package	1 piece	*Does not include the Y-Valve 2730 as shown in the picture.	
Part Number	V-861721		

The following parts are suitable for these products:

Vyntus™ ECG



All Accessories

<b>Disposable ECG electrodes, BlueSensor R, for stress ECG</b>		For Device	
		Vyntus ECG CardioSoft ECG	
Package	25 pieces	•	
Package	31439-001		
<b>Accu NiMH 2850 mAh, AA, 1.2 V</b>		For Device	
		Vyntus ECG	
Package	1 piece	*Accu ≥ 2500 mAh	
Package	V-806520		
<b>Carry bag Vyntus™ ECG</b>		For Device	
		Vyntus ECG	
Package	25 pieces		
Package	V-841234		

The following parts are suitable for these products:

CardioSoft EKG CAM-USB

All Accessories

<b>Disposable ECG electrodes, BlueSensor R, for stress ECG</b>		For Device	
		Vyntus ECG CardioSoft ECG	
Package	25 pieces	•	
Part Number	31439-001		
<b>Leadwire 10/set, grabber, AHA, for CAM-USB</b>		For Device	
		GE CardioSoft ECG: CAM-USB	
Package	1 set	*Set of 10 leadwires AHA with grabber; different lengths.	
Part Number	16822		

The following parts are suitable for these products:

Blood Pressure Accessories for VIAprint | Ergoselect Series | Suntech Tango

All Accessories

<b>Blood pressure cuff, Standard, 13 cm wide, with tube 130 cm, arm circumference 24 - 32 cm</b>		For Device	
		VIAprint 150 Ergoselect 600	
Package	1 piece	•	
Package	34589-001		
<b>Blood pressure cuff, Large, 15.5 cm wide, with tube 130 cm, arm circumference 32 - 42 cm</b>		For Device	
		VIAprint 150 Ergoselect 600	
Package	1 piece	•	
Package	30242-001		
<b>Blood pressure cuff, Small (children), 9 cm wide, with tube 130 cm, arm circumference 17 - 26 cm</b>		For Device	
		VIAprint 150 Ergoselect 600	
Package	1 piece	•	
Package	30243-001		
<b>Orbit-K BP cuff, Adult Small, Suntech, arm circumference 18 - 27 cm</b>		For Device	
		Tango M2	
Package	1 piece	•	
Package	25387-001		
		*For Tango and newer Lode ergometer with blood pressure units needing an ECG trigger.	





<b>Orbit-K BP cuff, Adult, Suntech, arm circumference 25 - 35 cm</b>		For Device	
		Tango M2	
Package	1 piece	•	
Part Number	25388-001		
		*For Tango and newer Lode ergometer with blood pressure units needing an ECG trigger.	
<b>Orbit-K BP cuff, Adult Plus, Suntech, arm circumference 27 - 40 cm</b>		For Device	
		Tango M2	
Package	1 piece	•	
Part Number	25389-001		
		*For Tango and newer Lode ergometer with blood pressure units needing an ECG trigger.	
<b>Orbit-K BP cuff, Adult Large, Suntech, arm circumference 32 - 44 cm</b>		For Device	
		Tango M2	
Package	1 piece	•	
Part Number	25390-001		
		*For Tango and newer Lode ergometer with blood pressure units needing an ECG trigger.	
<b>Orbit-K Microphone, 18" Suntech, for Tango M2 only</b>		For Device	
		Tango M2	
Package	1 piece	•	
Part Number	25386-001		






The following parts are suitable for these products:

Ergometer Accessories for VIAsprint 150P

**All Accessories**





<b>Pediatric saddle</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Package	26962-001		
<b>Racing saddle</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Package	26957-001		
<b>Quick release for saddle</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Package	26958-001		
<b>Horizontal saddle adjustment</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Package	26956-001		





<b>Adjustable pedal cranks</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Part Number	26959-001		
<b>Triathlon handlebar</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Part Number	26960-001		
<b>Anti tipping plate</b>		For Device	
		VIAsprint 150P	
Package	1 piece	•	
Part Number	26963-001		

The following parts are suitable for these products:

Vyntus™ WALK

All Accessories

<b>Adapter cable Nonin WO2 to DB9 Pin (3150 Oximeter)</b>		For Device	
		Vyntus WALK	
Package	1 piece	•	
Package	28439-001		
<b>SpO<sub>2</sub> Soft Sensor (3150 Oximeter), WO2 connector, 30 cm, size S</b>		For Device	
		Vyntus WALK	
Package	1 piece	•	
Package	22011-001		
<b>SpO<sub>2</sub> Soft Sensor (3150 Oximeter), WO2 connector, 30 cm, size M</b>		For Device	
		Vyntus WALK	
Package	1 piece	•	
Package	22010-001		
<b>SpO<sub>2</sub> Soft Sensor (3150 Oximeter), WO2 connector, 30 cm, size L</b>		For Device	
		Vyntus WALK	
Package	1 piece	•	
Package	22009-001		

<b>SpO<sub>2</sub> Fingerclip (3150 Oximeter), WO2 connector, 30cm</b>		For Device	
		Vyntus WALK	
Package	1 piece	•	
Part Number	29836-001		
<b>SpO<sub>2</sub> Flex Sensor Adult (3150 Oximeter), WO2 connector, 30 cm</b>		For Device	
		Vyntus WALK	
Package	1 piece	•	
Part Number	22008-001		
<b>Earclip SpO<sub>2</sub> sensor with 1 m cable</b>		For Device	
		Vyntus WALK	
Package	1 piece	• Requires additionally V-440044 Adapter cable Nonin WO2 to DB9 Pin (3150 Oximeter)	
Part Number	28896-001		
<b>Wraps for Adult Flex SpO<sub>2</sub> sensor, 8000JFW, disposable</b>		For Device	
		Vyntus WALK	
Package	25 pieces	•	
Part Number	22012-001		

Cardiopulmonary Exercise Testing (CPET),  
ECG and 6-Minute-Walk-Test

All Accessories

<p><b>Wrist band reusable 8-28 cm (3150 oximeter)</b></p>		For Device	
		Vyntus WALK	
Package	1 piece	<p>.</p>	
Package	22989-001		

# Index


Adapter .....	15	Pneumotach.....	13
Adapter .....	15	Protective cover .....	15
Adapter .....	20	Saddle .....	40
Adapter .....	20	Saddle.....	40
Adapter .....	25	Saddle.....	40
Adapter .....	32	Saddle.....	40
Adapter .....	34	Sealing .....	13
Adapter cable .....	42	Sealing .....	17
Anti tipping plate.....	41	SpO <sub>2</sub> extension cable.....	28
Bag.....	35	SpO <sub>2</sub> Reflect Sensor Holder .....	29
Battery.....	36	SpO <sub>2</sub> Sensor.....	28
Blood pressure cuff.....	38	SpO <sub>2</sub> Sensor.....	28
Blood pressure cuff.....	38	SpO <sub>2</sub> Sensor.....	29
Blood pressure cuff.....	38	SpO <sub>2</sub> Sensor.....	29
Blood pressure cuff.....	38	SpO <sub>2</sub> Sensor.....	29
Blood pressure cuff.....	39	SpO <sub>2</sub> Sensor.....	30
Blood pressure cuff.....	39	SpO <sub>2</sub> Sensor.....	42
Blood pressure cuff.....	39	SpO <sub>2</sub> Sensor.....	42
Calibration syringe.....	20	SpO <sub>2</sub> Sensor.....	42
Calibration syringe.....	20	SpO <sub>2</sub> Sensor.....	43
Carry bag .....	36	SpO <sub>2</sub> Sensor.....	43
Clamp.....	17	SpO <sub>2</sub> Sensor.....	43
Clamp.....	17	Tube .....	14
Cranks.....	41	Tube .....	14
ECG electrodes.....	36	Tube .....	32
ECG electrodes.....	37	USS Module .....	14
Elbow .....	13	USS Module .....	14
Elbow .....	13	USS Module .....	15
Filling assembly .....	35	Valve .....	16
Filter .....	10	Valve .....	16
Filter .....	10	Valve .....	17
Filter .....	10	Valve .....	34
Filter .....	10	Volume sensor .....	32
FPV block .....	16	Volume sensor .....	33
FPV block .....	16	Volume sensor .....	33
Gas.....	19	Volume sensor .....	34
Gas.....	19	Wraps .....	28
Gas.....	19	Wraps .....	43
Gas.....	31	Wrist band.....	44
Gas.....	31		
Gas.....	31		
Handlebar .....	41		
Headgear .....	25		
Headgear .....	25		
Headgear .....	26		
Head Support.....	35		
Leadwires.....	37		
LED holder .....	32		
LED holder .....	34		
Mask.....	24		
Mask.....	24		
Mask.....	24		
Mask.....	24		
Mask.....	25		
Mask.....	26		
Mask.....	26		
Microphone .....	39		
Mouthpiece .....	11		
Mouthpiece .....	11		
Mouthpiece .....	12		
Mouthpiece .....	12		
Mouthpiece .....	26		
Nose clip .....	11		
Nose clip .....	11		
Nose clip .....	27		
Nose clip .....	27		
Oxygen sensor.....	18		
Oxygen sensor.....	33		
Pneumotach.....	12		
Pneumotach.....	12		

All illustrations, descriptions and technical specifications are subject to change without prior notice. Please note that all products, services, or features of products and services may not be available in your local area. Please check with your local Vyair representative.

\* Vyair Medical is not the legal manufacturer.

#### GLOBAL HEADQUARTERS

Vyair Medical, Inc.  
26125 N. Riverwoods Blvd.  
Mettawa, IL 60045  
USA

 Vyair Medical GmbH  
Leibnizstrasse 7  
97204 Hoechberg  
Germany  
+49 931 4972-0 tel

For US use only.

© 2023 Vyair Medical, Inc. or one of its affiliates. All rights reserved. Vyair, the Vyair Medical logo and all other trademarks are trademarks or registered trademarks of Vyair Medical, Inc. or one of its affiliates. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labeling.  
VYR-US-2300061 V1.0

