

# Choose Vyaire™ Medical solutions for your durable medical equipment needs.

## Vyaire Medical is breathing new life into respiratory care.

While our name may be new to you, our products have been helping you save lives and provide care for over 65 years. Today, we manufacture and market more than 27,000 unique products for the diagnosis, treatment and monitoring of respiratory conditions across the continuum of care.

Vyaire was born of innovative brands that include AirLife®, Vital Signs®, Bear Medical Systems, Pulmonetic Systems, Sensormedics™, Viasys™ and many others; and we bring decades of industry leadership, clinical experience and unrivaled expertise.

**Now, as a single global “breathing company,” we are focused on improving patient outcomes and increasing value for our customers.**

### LTV® 1150 ventilator



The LTV® ventilator series, with greater than 85,000 units sold globally, has been providing ventilator support for pediatric to adults for since 1999. The LTV® 1150 ventilator provides portable, advanced ventilation for adult and pediatric patients at home or a post-acute care facility. At only 14 lbs, the ventilator provides a wide range of ventilation therapies to meet demanding patient needs, including volume control, pressure control, pressure support and spontaneous breath types.

### HMEs




Lightweight heat and moisture exchangers (HMEs) and heat moisture exchanger with filter (HMEF) provide low resistance to flow and are treated to provide bacteriostatic and hygroscopic (HCH) benefits. Our portfolio also includes bypass heat and moisture exchangers as well as the Edith HME, which meets requirements for adult and pediatric patients in various critical care environments.

### Suction kits



Optimize suctioning while reducing trauma. Our Tri-Flo™ catheters and catheter kits offer high suctioning efficiency in many catheter sizes and configurations to fulfill both clinician and patient needs.

### Oxygen delivery



Our oxygen therapy disposables help meet nearly every oxygen therapy need across the care continuum. We offer an extensive range of AirLife® nasal cannulas, oxygen masks and accessories.

### Closed suction



Our closed suction catheters and Verso® adapters facilitate access to patient airways while helping to maintain continuous ventilation.

### Trach care kits



Sterile, single-use tracheostomy care and cleaning kits include complete trach cleaning trays and trach care kits. The trays include three rigid basins for cleaning, soaking and rinsing; pre-sewn trach dressing; and various other supplies for optimizing trach care.

### Modudose™ unit dose solutions



These easy-to-use unit dose solutions contain sterile water, normal saline or half-normal saline solution and are available in different sizes to accommodate your inhalation and irrigation needs.

## PureSom headgear and chinstrap



Keep a cool head with these CPAP accessories. PureSom headgear and chinstraps are made with Breathe-O-Prene®, a breathable material that allows air to flow, keeping your patients cool and comfortable. Our universal headgear can be used with our CPAP series and other vendor masks.

## Nox-T3 portable sleep monitor



Compact and lightweight home sleep monitoring. The NOX-T3 portable sleep monitor is an easy-to-use device for sleep diagnosis that includes the ability to monitor the frequency and intensity of bruxism and titration of oral appliances for the treatment of sleep apnea and snoring. Combined with Noxturnal software, the monitor lets you view and listen to the snoring signal, which helps optimize work flow and enhance overall care.

## Portable spirometry



Hold simple spirometry in the palm of your hands. Our hand-held spirometers range from screening devices for chronic obstructive pulmonary disease (COPD) to portable, pocket-sized spirometers, to e-Diary applications. Additionally, our applications help comply with ATS/ERS standards, offer pre- and post-bronchodilator measurements, and provide support with a choice of predicted values and with quality checks.

For more information, please contact us at **+1.833.327.3284**.

## VYAIRE MEDICAL

26125 N. Riverwoods Blvd.  
Mettawa, IL

[vyaire.com](http://vyaire.com)

**For U.S. distribution only.**

Trademarks are the property of their respective owners.

© 2017 Vyaire. Vyaire, the Vyaire Logo and all other trademarks are property of Vyaire Medical, Inc. | CF/805/17/0139 (0118)

## ApneaRX® sleep therapy device



Help your patients sleep better. The ApneaRX® sleep therapy device provides cost-effective therapy for sleep apnea. An easy-to-use, comfortable and effective boil and bite oral appliance, it can be easily fit to patients. It is titrated in 1 mm increments without any other tools.

## PAP masks



Our assortment of sleep masks combine technology and comfort for high performance and a restful night of sleep. Aura nasal mask is the world's first nasal mask featuring custom fit technology and an active headgear connector to provide comfort with a reliable seal. IQ nasal mask features an adjustable sealing surface that provides a tailored fit in just seconds. MiniMe 2™ closely fits and maximizes comfort for your smallest patients. And the Mojo™ full face mask features a contoured shell that follows facial curves and has a skin-soft gel cushion to help seal masks and minimize skin irritation.

## Peak flow meter and e-Diary



Take spirometry testing home. Our peak flow meter and e-Diary are innovative, compact spirometry devices that enable patients to conduct monitoring at home. The MicroPeak precision peak flow meter comes with the ATS or EU scale. Immersible in water for disinfecting, the flow meter does not need to be dismantled, which can compromise the internal mechanism. The Micro Diary can record PEF, FEV1, FVC or PEV6 and loop data with symptom score and questions.

**vyaire**<sup>™</sup>  
M E D I C A L