

# ApneaRx<sup>®</sup> titratable “transition” oral appliance

An easy, affordable sleep apnea  
treatment option



**A boil and bite oral appliance that adjusts in 1 mm increments and complies with the new Medicare guidelines for oral appliance therapy.**

The ApneaRx<sup>®</sup> oral appliance from CareFusion was specifically designed to provide cost-effective therapy for mild to moderate sleep apnea. A defining feature of this boil and bite oral appliance is that the provider or patient can easily fit and titrate in 1 mm increments without screws, plastic bands or other accessories. The ApneaRx appliance provides an additional treatment option for patients with mild to moderate sleep apnea.



## Benefits

- Sleep MDs receive a transition device that may help them determine patients who are responsive and non-responsive to custom oral appliance therapy.
- Sleep dentists can provide a quick temporary fix to a broken custom oral appliance. As a trial device, the appliance can help undecided patients determine whether to pursue custom oral appliance therapy.
- Durable medical equipment suppliers can provide a transition treatment option for failed continuous positive airway pressure (CPAP) therapy.
- Patients receive affordable, easy-to-use, comfortable and effective transition therapy for mild to moderate sleep apnea.
- Frequent travelers and patients with respiratory problems may have an alternative to CPAP.
- Patients requiring a cost-effective solution have a custom appliance option.
- Snoring patients without sleep apnea can use the ApneaRx oral appliance to help reduce their snoring.
- Technicians or patients can fit the appliance, and patients adjust it as needed.



## Features

- Advanced design uses a simple unaided mechanism to maintain a set jaw position.
- Comfort-Fit™ thermally molds to the patient's mouth in 10 minutes.
- Posi-Lock™ sets and locks the jaw position at any time.
- Micro-Fit™ provides easy-to-read 1 mm jaw adjustment to increase comfort and effectiveness.
- V-Flow™ assures full airflow and device alignment to enhance stability and fit.

For more information or to place an order, contact CareFusion Ordering and Customer Support at [800.582.7458](tel:800.582.7458) or [tiasales@carefusion.com](mailto:tiasales@carefusion.com).

CareFusion  
Yorba Linda, CA

[carefusion.com](http://carefusion.com)

© 2012 CareFusion Corporation or one of its subsidiaries. All rights reserved. CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its subsidiaries. ApneaRx is a trademark or registered trademark of Apnea Sciences Corporation or one of its subsidiaries. All other trademarks are property of their respective owners. RC760 (1012/1000)

