



GreenLight™
Laryngoscope handles and blades

Designed to make reprocessing easy

- Stainless steel and plastic ISO standard connection single-patient-use blades reduce time and effort associated with reprocessing of reusable blades
- Standard AA batteries* (*alkaline or lithium*) used in traditional handle for easy replacement (*single CR-123 lithium battery used in stubby handle*)
- With our fiber optic system there is no longer a need for bulb replacement

Designed for performance

- Ergonomic and over-molded handles are one of the lightest in the market
- Handle has a durable, slip-resistant grip and a wide, sturdy base
- Bright LED light and no back-flash glare
- Easy to remove battery and light cartridge for quick sterilization between cases. Handle is autoclavable

Designed to last

- Efficient technology extends battery life beyond that of a standard incandescent bulb light source
- Five-year warranty on all handles



Standard handle with single-patient-use stainless steel blade



Stubby handle with single-patient-use stainless steel blade

* Batteries should be disposed of as per local regulations

Part no.	Description	Handle type	Case qty
4558GSP	GreenLight™ ISO standard laryngoscope handle (Santoprene over-molded handle) - reusable	Standard	1
4559GSP	GreenLight™ ISO standard laryngoscope handle (Santoprene over-molded handle) - reusable	Stubby	1
4559LED	Replacement battery/LED cartridge for GreenLight™ (ISO green standard laryngoscope handle)	Standard	1
4559LEDS	Replacement battery/LED cartridge for GreenLight™ (ISO green standard laryngoscope handle)	Stubby	1



Standard handle



Stubby handle



Stainless steel blades



Plastic blades

Part no.	Description	Blade type	Blade size	Case qty
4602	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Mac	2	20
4603EU	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Mac	3	20
4604EU	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Mac	4	20
4610	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	0	20
4611EU	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	1	20
4612	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	2	20
4613	GreenLight™ ISO standard, stainless steel, SPU laryngoscope blade	Miller	3	20
4002GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Mac	2	20
4003GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Mac	3	20
4004GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Mac	4	20
4010GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	0	20
4011GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	1	20
4012GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	2	20
4013GSP	GreenLight™ ISO standard, plastic, SPU laryngoscope blade	Miller	3	20



For more information, please contact your local Vyairē representative.

GLOBAL HEADQUARTERS

Vyairē Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Vyairē Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



AUSTRALIAN SPONSOR

Vyairē Medical Pty Ltd
Level 5, 7 Eden Park Drive
Macquarie Park, NSW, 2113
Australia

vyaire.com

For Global Distribution.

© 2019 Vyairē. Vyairē and the Vyairē logo and all other trademarks or registered trademarks are property of Vyairē Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-1900115