



AirLife™ volumetric incentive spirometer

Deep breathing therapy is simplified with AirLife™ volumetric incentive spirometers. Their intuitive design prompts patients to correctly perform and monitor their own breathing exercises, even without direct supervision. Both patients and clinicians will appreciate:

Convenient handle

Place right on the bed rail for safe storage and easy accessibility.

Enhanced legibility

Dual-sided graduations allow for fast and easy interpretation of results.

Durable, compact design

The sturdy plastic resists accidental breakage for improved patient safety and reduced cost to hospitals.

The AirLife volumetric incentive spirometer portfolio offers these quality features

Highly legible calibrations

Dual-sided calibrations clearly identify achieved volumes of up to 4,000 mL capacity. Clinicians and patients can easily read results from either side of the unit.

Patient goal indicator

Easily adjustable indicator provides a volumetric goal for adequate lung exercise and allows patients to help monitor their own progress.

Easy to grip handle

Handle is easy to grip and helps ensure that patients use the spirometer in an upright position for proper therapy. Handle can also be placed onto bed rail for safe storage and easy accessibility.



Flexible tubing/swivel mouthpiece

Mouthpiece allows the patient to easily adjust for use in various positions. Mouthpiece holder provides a place to store the mouthpiece when not in use.

Inspiratory flow indicator

Indicator encourages patients to maintain a slow, proper rate of inspiration. Slow inspiratory flow rate promotes even distribution of inspired volume for enhanced therapy.

One-way valve option

Provides a reminder for the patient to inhale rather than exhale, ensuring correct technique, even when you can't be at the bedside.

Choose from a 2,500 or 4,000 mL adult size or the 2,500 mL pediatric version with colorful, eye-catching graphics.

Ordering information

Cat. no.	Description
4,000 mL capacity	
001901A	With one-way valve
001902A	Without one-way valve
2,500 mL capacity	
001903A	With one-way valve
001904A	Without one-way valve
2,500 mL capacity (pediatric version)	
001905A	With colorful targets and one-way valve

All spirometers are packaged 12 per case.

For more information, contact your local Vyair sales representative today, or visit our website.

GLOBAL HEADQUARTERS

Vyair Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



Vyair Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



AUSTRALIAN SPONSOR

Vyair Medical Pty Ltd
Suite 5.03, Building C
11 Talavera Road
Macquarie Park, NSW 2113
Australia

vyaire.com

REGONLY For global distribution.

© 2020 Vyair. Vyair and the Vyair logo are trademarks or registered trademarks of Vyair Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-2000412

vyaire[™]
M E D I C A L