

vyaire™
MEDICAL



Misty Fast™
Disposable small-volume nebuliser



*“Support your patient’s
next breath while
reducing treatment time”*

Misty Fast™ disposable small-volume nebuliser (SVN)

The Misty Fast™ (SVN) is a small volume, high-efficiency air entrainment nebuliser that delivers comparable results to the Misty Max 10™ nebuliser, but in half the time. The one-piece jet design provides efficient, high output rates, for fast treatment times and consistent medication delivery. Eight configurations accommodate your patient’s aerosolisation needs, and all include our crush-resistant O₂ tubing.



The Misty Fast™ SVN offers the following key benefits:

- Treatment requires half the time of the AirLife Misty Max 10 SVN
- Proprietary swivel mouthpiece helps enhance patient positioning
- Unique, integrated jet/baffle design provides medication delivery consistency and reliability
- Clear base with visible graduations facilitates medication loading and aerosol delivery visualisation

Misty Fast SVN (50 per case)

Cat no.	Product description
FN2431	Misty Fast SVN and 2.1 m O ₂ tubing
FN2433	Misty Fast SVN with adult aerosol mask and 2.1 m O ₂ tubing
FN2434	Misty Fast SVN with baffled tee swivel mouthpiece and 2.1 m O ₂ tubing
FN2438	Misty Fast SVN with baffled tee swivel mouthpiece, 15.2 cm flextube and 2.1 m O ₂ tubing
FN2439	Misty Fast SVN with baffled tee swivel mouthpiece, 15.2 cm flextube and 4.3 m U/Connect-It® O ₂ tubing
FN2444	Misty Fast SVN with pediatric aerosol mask and 2.1 m O ₂ tubing
FN2445	Misty Fast SVN with baffled tee swivel mouthpiece and 2.1 m U/Connect-It O ₂ tubing
FN2446	Misty Fast SVN with baffled tee swivel mouthpiece, 15.2 cm flextube and 2.1 m U/Connect-It O ₂ tubing



For more information, please contact your local Vyairé representative.

GLOBAL HEADQUARTERS

Vyairé Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Vyairé Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA



AUSTRALIAN SPONSOR

Vyairé Medical Pty Ltd
Level 5, 7 Eden Park Drive
Macquarie Park, NSW, 2113
Australia

vyaire.com

For Europe, Asia, Australia/New Zealand, Middle East and Africa distribution only.

© 2019 Vyairé Medical, Inc. or one of its affiliates. All rights reserved. Misty Fast, Vyairé and the Vyairé Medical logo are trademarks or registered trademarks of Vyairé Medical, Inc or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-INT-1900142