

# New & Improved ApexPro™ FH adapter

We heard you! And as we improved our adapter design, we took all of your feedback into consideration, making it smaller, easier to use, and more functional.



The new ApexPro FH Adapter features:

- Smaller design making it easier for patients to move about the hospital while being monitored walking with the device.
- Closure plug keeps the ports clean while using a 5 lead and helps lengthen life of device by protecting ports.
- Tested for chemical resistance of cleaning and disinfection agents instructed in product IFU.
- Compatibility tested within the hospital environment for ApexPro FH system signal accuracy according to AAMI 60601-2-27".



## GLOBAL HEADQUARTERS

Vyair Medical, Inc.  
26125 North Riverwoods Blvd  
Mettawa, IL 60045, USA

[vyaire.com](http://vyaire.com)

**Rx ONLY** For U.S. distribution only.

©2019 Vyair. Vyair and the Vyair Logo are trademarks of Vyair Medical, Inc. | VYR-US-1900083

**vyaire**<sup>™</sup>  
M E D I C A L