

Fenom Pro[®] Asthma Monitor

Technical Data Sheet

Fenom Pro[®]



Dimensions and weight	Height: 145 mm, Width: 230 mm, Depth: 140 mm Weight (including handpiece): 2.4 kg (5.3 lbs)
Electrical data	Device power consumption: < 20 VA Power supply mains voltage: 100–240 V ~50-60 Hz
Operating Range	Ambient temperature between 15-30°C (59-86°F) Relative humidity between 20-80% Ambient pressure between 80-106 kPa
Measurement range	10 ppb to 200 ppb
Measurement time	Results will display in 28 seconds
Precision	< 5 ppb or 10% of concentration
Accuracy	±5 ppb or ±10% of concentration
Limit of detection	10 ppb
Exhalation parameters	Exhalation time: 10 seconds Exhalation pressure is between 15-25 cm-water Exhalation flow rate is between 45-55 mL/s
Rechargeable battery	Capacity: > 15 tests over 6 hours on a fully charged battery Charging time: 4 hours
Device life and service	Fenom Pro device life is 5 years. Periodic service by your local representative is required.
Manufacturer information	Spirosure, Inc. 7020 Koll Center Parkway, Suite 110, Pleasanton, CA 94566-3107

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA

Rx ONLY For U.S. distribution only.

Trademarks are the property of their respective owners.

© 2020 Vyaire. Vyaire and the Vyaire logo are property of Vyaire Medical, Inc., or one of its affiliates. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-US-2000006