


OxyJet®

Provides supplemental oxygen for patients in need during use of the vPEP® or PocketPEP®

Designed by clinicians for clinicians and their patients

The OxyJet® is an innovative venturi adapter to provide supplemental oxygen during use of the vPEP®. This device must be used with vPEP® or PocketPEP® (with additional T-piece adapter, item 557V).



vPEP® combined with OxyJet®

OxyJet features and benefits*:

- ✓ Compact and ergonomic design
- ✓ Provides adjunct positive airway pressure
- ✓ Helps decrease work of breathing
- ✓ Easy-to-clean** to reduce risk of infection
- ✓ Easy-to-read manometer for monitoring pressure so oxygen flow rate can be adjusted as needed
 - Comes with standard 7' tubing

* When used with the vPEP or PocketPEP device

** The OxyJet version meets Cystic Fibrosis Foundation recommended cleaning guidelines¹

OxyJet® specifications:

- Clean with soapy water, dishwasher, boiling or microwave steam bags
- PVC and DEHP free
- Not manufactured from latex natural rubber
- Single patient, multiple use device



Item name	Item number	Quantity per case
OxyJet with 7 ft tubing	561V	1 case of 10 each

For more information, please contact your local Vyaire representative or visit vyaire.com.

REFERENCES

1. Cystic Fibrosis Foundation. Nebulizer Care at Home. Cystic Fibrosis Foundation. <https://www.cff.org/Life-With-CF/Treatments-and-Therapies/Medications/Nebulizer-Care-at-Home/>

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045, USA

vyaire.com

Rx ONLY For U.S. distribution only.

OxyJet, vPEP and PocketPEP are registered trademarks of D.R. Burton Healthcare.
© 2019 Vyaire. Vyaire and the Vyaire logo are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. VYR-US-1900240

vyaire[™]
M E D I C A L