



AVEA™ ventilator
Disposable Expiratory Filter / Watertrap

The HEPA filtration element filters the expired gas to 0.3 microns. With a 99.999+% bacterial and viral efficiency rating, and the heated environment, the expiratory flow sensor minimizes the risk of contamination, while providing an exceptional environment for your staff and patients.

Resistance comparison		
Filter type	Ventilator	Published resistance
Disposable Expiratory Filter / Watertrap (Part number 11790)	AVEA ventilator	< 1.0 cm H ₂ O / 60 LPM Flow
Pall with reusable expiratory filter	AVEA ventilator	2.8 cm H ₂ O / 100 LPM Flow 1.75 cm H ₂ O / 60 LPM Flow

Description	Part number	Quantity
AVEA ventilator single patient use expiratory filter / water trap	11790	Case of 12

Benefits of using the AVEA Disposable Expiratory Filter / Watertrap:

- Single patient use
- No component assembly required
- Easy condensation management
- 60% lower resistance than the reusable configuration

⚠ WARNING: U.S. Federal Law restricts this device to sale by or on the order of a physician.

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA

vyaire.com

Rx ONLY For U.S. distribution only.

© 2020 Vyaire. Vyaire, the Vyaire Logo and all other trademarks are property of Vyaire Medical, Inc. VYR-US-2000073

vyaire[™]
M E D I C A L