


Vyaire X2 Solutions


Spacelabs telemetry adapter and direct connect

The new adapters and leadwires expands from Spacelabs bedside monitors to the telemetry system.

- Adapter allows compatibility to full 5-lead Vyaire X2 single-patient-use leadwire family
- Small and lightweight design maximizes patient experience
- Direct connect leadwires connect directly to the telemetry transmitter if there is no need to move between telemetry and bedside monitoring
- Enables use of Multi-Link™ reusable leadsets
- Tested for chemical resistance of cleaning and disinfection agents listed in product IFU.
- Compatibility tested within the hospital environment for AriaTele™ system signal accuracy according to AAMI 60601-2-27



2090601-712



2090601-070

SKU	Item description	Monitor platform	Product type
2090601-070	Spacelabs Aria 5-lead Adapter	Spacelabs AriaTele	Adapter
2090601-712	Spacelabs Aria 5-lead direct connect SPUL 60cm – grabber	Spacelabs AriaTele	Leadwire
2090601-722	Spacelabs Aria 5-lead direct connect SPUL 60cm- snap	Spacelabs AriaTele	Leadwire

GLOBAL HEADQUARTERS

Vyaire Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045, USA



Vyaire Medical Oy
Kuortaneenkatu 2
FI-00510 Helsinki, Finland



AUSTRALIAN SPONSOR

Vyaire Medical Pty Ltd
Suite 5.03, Building C
11 Talavera Road
Macquarie Park, NSW 2113
Australia

vyaire.com

Trademarks are the property of their respective owners.

RE-ONLY For global distribution only

© 2020 Vyaire. Vyaire and the Vyaire logo are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-GBL-2000537

vyaire[™]
M E D I C A L